


DR400C Patient Scale

Accurate and reliable, Detecto's economical DR400C platform scale is lightweight and portable, perfect for mobile clinics and home care nurses. The remote indicator features a large 1" display, units conversion, and tare. Battery powered so you can have a scale anywhere you need it.



SPECIFICATIONS

- Capacity:** 400 lb x .5 lb / 181 kg x .2 kg
- Display:** 4 digit, 7 segment, liquid crystal display, 1" / 25 mm high
- Keypad:** 4 keys: lb/kg, Tare/Hold, On, Off
- Platform:** 12" square, painted steel, non-skid surface
- Units:** Pounds, kilograms
- Load Cells:** 1
- Power:** 9 volts DC supplied by 6 "AA" alkaline batteries or 115 volt AC adapter (included)

OPTIONAL CARRYING CASE AVAILABLE

Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.