



UNIVERSAL GLOVE BOX HOLDERS

MADE IN THE USA

Detecto Lifetime Guarantee

**STAINLESS OR
MILD STEEL**

**3 SIZES
AVAILABLE**



Detecto's model GH2SS (pictured above) holds two glove boxes.

DETECTO'S UNIVERSAL GLOVE BOX HOLDERS

Detecto's model GH Glove Box Holders are available in three different sizes and your choice of stainless steel or mild steel with white powder coat finish. Fit one glove box in the small GH1/GH1SS, two boxes side by side in the GH2/GH2SS, and three boxes side by side in the GH3/GH3SS. The glove box holders carry a Detecto lifetime guarantee and are quality made in the USA.

The units install in seconds via screwholes and the GH1/GH1SS also features a snap-on clip in back. The open design allows you to see what size or type of glove or tissue is being dispensed.



DETECTO LIFETIME GUARANTEE



The GH's are universal enough to be used for tissue box holders where counter space is conservative.



The model GH1SS holds one box snugly and accommodates most glove box types.

Mild Steel Glove Box Holders

MODEL	DIMENSIONS	SHIPPING WEIGHT
GH1 (Holds 1 Box)	10.5" H x 6.5" W x 4.25" D / 27 cm H x 17 cm W x 11 cm D	3 lb / 1.4 kg
GH2 (Holds 2 Boxes)	10" H x 12" W x 4" D / 25 cm H x 30 cm W x 10 cm D	3 lb / 1.4 kg
GH3 (Holds 3 Boxes)	10" H x 18" W x 4" D / 25 cm H x 46 cm W x 10 cm D	3 lb / 1.4 kg



Stainless Steel Glove Box Holders

MODEL	DIMENSIONS	SHIPPING WEIGHT
GH1SS (Holds 1 Box)	10.5" H x 6.5" W x 4.25" D / 27 cm H x 17 cm W x 11 cm D	3 lb / 1.4 kg
GH2SS (Holds 2 Boxes)	10" H x 12" W x 4" D / 25 cm H x 30 cm W x 10 cm D	3 lb / 1.4 kg
GH3SS (Holds 3 Boxes)	10" H x 18" W x 4" D / 25 cm H x 46 cm W x 10 cm D	3 lb / 1.4 kg



Cardinal Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By:

