



# Digital Beam Scale

## Model DS5100

Step up to a faster and easier way to weigh patients with the Model DS5100 Digital Beam Scale. The ideal scale for a busy physician's office, the DS5100 displays the patient's weight in just three seconds after they step on the large 15" x 13" low profile platform. Sliding counterweights on a beam and waiting for it to settle is a thing of the past.

There is no need to press the power button to turn on the scale, thanks to the Auto Start-up feature. Stepping on the platform turns on the DS5100, and the scale automatically shuts off when not in use.

The sturdy aluminum column houses the EZ-Slide height measurement rod. Enter the height using the scroll buttons, press the BMI button and the patient's Body Mass Index is instantly displayed in the right LCD.

**MAXIMUM CAPACITY**  
500 x 0.1 lb (200 x 0.05 kg)

**BMI CALCULATOR**  
Input data in FT/IN or CM

**ACCURACY**  
+/- 0.1 lb (0.05 kg)

**POWER**  
Four C Batteries (Included)

**PLATFORM SIZE**  
15" W x 13" D x 2" H

**BATTERY LIFE**  
Four Years or 50,000 Weighments

**DISPLAY**  
Dual 1" High LCD

**CONSTRUCTION**  
ABS and Metal Construction

**CLOCK**  
12 or 24 Hour Display

**WARRANTY**  
Three Years



*Featuring*  
**Motion Lock**  
*Technology*



The Model DS5100 features large, easy to read controls. Data is displayed in dual 1" high LCD's.