



DS5200 WAIST HIGH DIGITAL PHYSICIAN SCALE

Revolutionize your clinic  
with our cost-efficient  
waist-high solution.



# Waist High Digital Physician Scale

## KEY FEATURES

For today's busy physician offices, Doran's DS5200 Waist High Digital Physician Scale provides a faster, easier way to weigh and measure patients accurately — *at a price that cannot be beat!* Featuring a large slip-resistant platform, the DS5200 also has the lowest step up height in the market. Our *Touchless Start* feature allows you to acquire a patient's weight in as little as two seconds — without any buttons to press. Once the patient steps onto the platform, an accurate weight up to 600 pounds will be held on the display. When the patient steps off, the DS5200 is immediately ready for another weighing. With our *Units Lock* feature, your scale will never display an incorrect unit of measure.

Doran's exclusive *Weight Lock* feature ensures an accurate weight is displayed regardless of patient movement. Our *Reweigh* functionality also allows a second weight to be acquired while the patient remains on the platform.

The industry's only internal rechargeable battery is designed to last — providing at least 5500 weighments between charges. When a recharge is needed, just connect the USB port to any wall mount charger and be fully recharged in six hours.

Easily calculate *Body Mass Index* using the scroll buttons to enter the height, then press the BMI button. The patient's index is instantly displayed. The DS5200 is EMR ready and data can be transmitted via the included USB port using our communications option.

## DETAILS THAT MATTER

The DS5200's platform features the industry's lowest step-up height and durable wheels for enhanced patient comfort with easy mobility.



## SPECIFICATIONS

<b>Warranty:</b>	Three Years
<b>Power:</b>	Internal rechargeable battery. USB charging cable included. Recharges with standard 5V, 1A USB adapter (not included).
<b>Battery Life:</b>	5500 Weighments
<b>Maximum Capacity:</b>	600 x 0.2 lb (270 x 0.1kg)
<b>Dimensions:</b>	Platform Dimensions 14.0"W x 18.0"L x 1.35"H Overall Dimensions 14.0"W x 20.5"L x 35.5"H
<b>Display:</b>	1.25" High LED Display
<b>Displayed Units:</b>	lb, kg, Units Lock Feature
<b>Controls:</b>	Zero, Reweigh, Units, BMI, Up / Down Arrows
<b>BMI Calculator:</b>	Input data in FT/IN or CM
<b>Communications:</b>	EMR ready (with USB communications options)



## SMART DESIGN

The DS5200's elegant and ergonomically designed indicator is easy to read and operate. It offers robust functionality that makes fast and accurate weighments a snap.



DS5250 EYE LEVEL DIGITAL PHYSICIAN SCALE

Precision and Affordability.  
The ideal physician scale  
for weighing excellence.



# Eye Level Digital Physician Scale

## KEY FEATURES

For today's busy physician offices, Doran's DS5250 Eye Level Digital Physician Scale provides a faster, easier way to weigh and measure patients accurately — *at a price that cannot be beat!* Featuring a large slip-resistant platform, the DS5250 also has the lowest step up height in the market.

Our *Touchless Start* feature allows you to acquire a patient's weight in as little as two seconds — without any buttons to press. Once the patient steps onto the platform, an accurate weight up to 600 pounds will be held on the display. When the patient steps off, the DS5250 is immediately ready for another weighing. With our *Units Lock* feature, your scale will never display an incorrect unit of measure.

Doran's exclusive *Weight Lock* feature ensures an accurate weight is displayed regardless of patient movement. Our *Reweigh* functionality also allows a second weight to be acquired while the patient remains on the platform.

The industry's only internal rechargeable battery is designed to last — providing at least 5500 weighments between charges. When a recharge is needed, just connect the USB port to any wall mount charger and be fully recharged in six hours.

The sturdy and handsome aluminium column elegantly houses the easy-to-use height measurement rod. Simply extend the measurement arm and pull upwards for adult heights or lower it to measure a child. Easily calculate *Body Mass Index* using the scroll buttons to enter the height, then press the BMI button. The patient's index is instantly displayed. The DS5250 is EMR ready and data can be transmitted via the included USB port using our communications option.

## DETAILS THAT MATTER

The DS5250 is equipped with durable wheels for reliable mobility. Our optional arms package adds extra safety and makes portability even easier.



## SPECIFICATIONS

<b>Warranty:</b>	Three Years
<b>Power:</b>	Internal rechargeable battery. USB charging cable included. Recharges with standard 5V, 1A USB adapter (not included).
<b>Battery Life:</b>	5500 Weighments
<b>Maximum Capacity:</b>	600 x 0.2 lb (270 x 0.1kg)
<b>Dimensions:</b>	Platform Dimensions 14.0"W x 18.0"L x 1.35"H Overall Dimensions 14.0"W x 20.5"L x 49.75"H Overall Dimensions with Arms 23.8"W x 20.5"L x 49.75"H
<b>Display:</b>	1.25" High LED Display
<b>Displayed Units:</b>	lb, kg, Units Lock Feature
<b>Controls:</b>	Zero, Reweigh, Units, BMI, Up / Down Arrows
<b>BMI Calculator:</b>	Input data in FT/IN or CM
<b>Communications:</b>	EMR ready (with USB communications options)
<b>Height Rod Adjustment:</b>	Minimum 2' 0" (61 cm) Maximum 6' 6.7" (200 cm)

## SMART DESIGN

The DS5250's elegant and ergonomically designed indicator is easy to read and operate — with robust functionality that makes fast, accurate weighments a snap.

