

CHAD  
**drive**

100% US QUALITY  
INSPECTED!

**NEW!**



**CH4808-L-BLUE &  
CH4815-L-BLUE**

**870 Regulator**

For use with: M2 through  
E sized cylinders

**Available** **CHAD® Regulators**

**1 Pack**

**870 REGULATORS**

CH4808-L-BLUE 1/ea 0-8 lpm, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8

CH4815-L-BLUE 1/ea 0-8 lpm, .5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15

**FEATURES:**

- Brass components in high-pressure areas and relief valve zones.
- 20-micron, sintered-bronze inlet filter.
- FDA-recommended viton and brass seal washers.
- All flows tested to accuracy within +/- 15% tolerances.
- Meets ASTM-G175-03 promoted ignition tests.
- Ergonomic t-handle for easy tightening.
- Inlet pressure: 200 - 3,000 PSI.
- Outlet pressure: 50 PSI.
- **100% QC tested in the USA.**
- Two Year Warranty.



**18301G**

**870 Regulator**

For use with: M2 through  
E sized cylinders

**Available**  
**1 Pack**

**Regulators**

**870 REGULATORS**

18300G Barb Outlet, Pediatric, Green, 1/bx,  
Liter Flow Increments: .03, .06, .12, .25, .5,  
.75, 1, 1.5, 2, 2.5, 3

18301G Barb Outlet, Adult, Green, 1/bx  
Liter Flow Increments: .12, .5, 1, 1.5, 2, 2.5, 3, 4, 6, 8

18302G Barb Outlet, Adult, Green, 1/bx  
Liter Flow Increments: .5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15



**18301GM**

**870 Mini Regulator**

For use with: M2 through  
E sized cylinders

**Available**  
**1 Pack**

**870 MINI REGULATORS**

18301GM Barb Outlet, Adult, Green,  
(Has replaced 18301, 18311), 1/bx  
Liter Flow Increments: .12, .5, 1, 1.5, 2, 2.5, 3, 4, 6, 8

18302GM Barb Outlet, Adult, Green, 1/bx  
Liter Flow Increments: .5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15

**FEATURES:**

- New click style flow control.
- Lightweight uni-body design.
- Meets or exceeds accuracy standards for ASTM, American National and CGA.
- UL Approved.