

IV POLES

Standard and Deluxe Models



WE CARE LIKE FAMILY™

IV POLE

- Powder-coated steel construction
- Removable 2-hook hanger
- Sturdy 4-leg base with two locking casters
- Tool-free, adjustable height, chrome extension pole with locking collar: 46"– 81.5"
- Easy to assemble and store
- 4-hook hanger option also available (10554-4)

DELUXE IV POLE

- Chrome-plated steel pole
- Removable 2-hook hanger
- Tool-free adjustable height with locking collar
- 5-leg steel base
- 4-hook hanger option also available (10555-4hook)



**IV POLE
(2 HOOKS)**



**IV POLE
(4 HOOKS)**



**DELUXE IV POLE
(2 HOOKS)**



**DELUXE IV POLE
(4 HOOKS)**

Reorder No.	Description	Qty	HCPCS
10554-2hook	IV Pole - Economy 4 Legs with 2 Hooks	1/Cs	E0776
10554-4hook	IV Pole - Economy 4 Legs with 4 Hooks	1/Cs	E0776
10555-2hook	IV Pole - Deluxe 5 Legs with 2 Hooks	2 Bx	E0776
10555-4hook	IV Pole - Deluxe 5 Legs with 4 Hooks	2 Bx	E0776



Find Us On:



dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR

Call us at: (888) 396-2739 • www.dynarex.com