



WE CARE LIKE FAMILY™

# TOUCH SCREEN VITAL SIGNS PATIENT MONITOR with STAND



- Accurately monitors and records NIBP, pulse rate, SpO<sub>2</sub>, and temperature
- Stores up to 960 hours of trend data
- 7" anti-glare touchscreen with large font display provides improved visibility
- Interface features interactive pop-up menus and dialogue boxes
- Includes 2 cuff sizes: Adult (9.8"-13.8") and Adult Large (13"-18.5"), Masimo Pulse Oximeter, Covidien Temperature Sensor, and 100 Probe Covers
- Covidien sensor increases the accuracy of readings
- Comes with a height-adjustable rolling stand for increased maneuverability and portability
- Includes rechargeable lithium-ion battery
- Suitable for adult, pediatric, and neonatal patients when used with the appropriate sensors



SpO<sub>2</sub> Sensor



Temperature Sensor



Blood Pressure Cuff

Reorder No.	Description	Qty./Bx
10841-w/Stand	Vital Signs Patient Monitor, Touch Screen with Stand	1

www.dynarex.com • 888-396-2739  
Dynarex Corporation • 10 Glenshaw Street • Orangeburg, NY 10962

R210224



© 2021 Dynarex Corporation