



Ultrasound Gel



Medical Imaging Assistance

Ultrasound Gel

- Electroconductive gel media used to facilitate the ultrasonic pulsed echo imaging system
- Acoustically correct for a broad range of frequencies
- Aqueous, will not stain clothing or damage transducers
- Can be used with any equipment that requires an aqueous, conductive gel (for external use only)
- Fragrance free
- Bacteriostatic and non-irritating
- Non-greasy
- Easy refill



Reorder No.	Ultrasound Gel	Quantity
1241	8 oz Blue	12/cs
1243	1.3 gal (5 liter) Jug, Blue	4/cs
1245	8 oz Clear	12/cs
1247	1.3 gal (5 liter) Jug, Clear	4/cs

Already mobile?
Scan your smart phone here for more information.



Visit our website to see our DynaLube (#1250-1254) and our Stethoscopes (#7110-7111)



For more information contact your Dynarex sales representative
888-396-2739

*Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.