



OXYGEN CONCENTRATOR

The Dynarex Resp-O₂™ Oxygen Concentrator is designed to provide supplemental oxygen to help relieve breathing difficulties in respiratory patients. The 5-liter Oxygen Concentrator includes audible alarms, oxygen purity indicator (OPI), and user-friendly flow-meter knob to provide prescribed liters-per-minute. The 37-pound concentrator is easily transportable with accessible top handle and four rolling casters.

5 Liter

REORDER NO. 33950

PACKAGE QUANTITY
1/Case

- Features a flow range of 0.5-5 LPM
- Audible and visual alarms include Low Flow-Rate, Low/High Pressure, and Power Failure
- Manual reset button helps ensure uninterrupted oxygen delivery
- Designed for quiet operation with low noise output
- Easy-to-read control panel for optimal visibility

ALSO AVAILABLE:

Dry Bubble Humidifier – 400 ML • 6 PSI (HOME USE)

REORDER NO. 34221

Oxygen Humidifier Elbow Connectors with Tubing

REORDER NO. 34505

Oxygen Humidifier Straight Connectors with Tubing

REORDER NO. 34506

Sterile Water for Inhalation, USP – PREFILLED BUBBLE HUMIDIFIER

REORDER NO. 34300 – SIZE 350 ML

REORDER NO. 34301 – SIZE 550 ML



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



R210823

©2021 Dynarex Corporation



OXYGEN CONCENTRATOR ACCESSORIES

Sterile Water for Inhalation

Dynarex Resp-O₂™ Sterile Water supplies are disposable, single-use respiratory therapy bottles filled with water designed for patient inhalation. The water is sterilized to ensure patient safety and to help avoid cross-contamination.

- Reservoir bottles feature built-in oxygen tubing connectors
- Sterile, transparent adaptor body provides optimal visualization
- Audible alarm alerts clinician to air-flow constrictions

Bubble Humidifier

The Dynarex Resp-O₂™ Disposable Bubble Humidifier is designed to humidify supplemental oxygen for patients receiving oxygen therapy at home. Ideal for use with all types of oxygen concentrators with either face masks or cannulas, the Bubble Humidifier is intended for oxygen low flow rates up to 6 PSI. The transparent humidifier features clearly defined minimum and maximum lines for easy monitoring.

- Features pressure-relief valve with audible alarm for optimal safety
- Includes PVC diffuser to help minimize noise
- Disposable, single-patient use for optimal hygiene
- Not made with natural rubber latex or DEHP

Oxygen Humidifier Connectors

Dynarex Resp-O₂™ Oxygen Humidifier Connectors are designed to attach dry humidifier bottles or pre-filled bubble humidifiers to continuous-flow oxygen concentrators. They come with either sturdy oxygen tubing on one end and a male threaded grip elbow on the other end, or with star lumen oxygen tubing on one end and molded threaded straight connector on the other end. Connectors are universal and fit most humidifier bottles with a DISS connection.

- Transparent oxygen tubing allows for easy visualization
- Universal fit connects with most humidifier bottles and concentrators
- Ribbed connector provides a secure fit without damaging connection points
- Designed to attach to barb-style concentrator outlets
- Not made with natural rubber latex or DEHP



Sterile Water for Inhalation, USP

PREFILLED BUBBLE HUMIDIFIER

REORDER NO. 34300

SIZE
350 mL

PACKAGE QUANTITY
20/Case (with adaptor)

Sterile Water for Inhalation, USP

PREFILLED BUBBLE HUMIDIFIER

REORDER NO. 34301

SIZE
550 mL

PACKAGE QUANTITY
12/Case (with adaptor)

Dry Bubble Humidifier

400 ML • 6 PSI (HOME USE)

REORDER NO. 34221

PACKAGE QUANTITY
50/Case

Oxygen Humidifier Elbow Connectors with Tubing

REORDER NO. 34505

PACKAGE QUANTITY
10/Case

Oxygen Humidifier Straight Connectors with Tubing

REORDER NO. 34506

PACKAGE QUANTITY
50/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.

©2021 Dynarex Corporation



R210823