


Products for the Laboratory Industry



Shoe Covers

Reorder No.	Description	Quantity
2131	Shoe Cover - Non-Conductive, blue	150 pr/Cs
2132	Shoe Cover - Non-Conductive, Non-Skid, blue	150 pr/Cs
2134	Shoe Cover - Non-Conductive, Non-Skid, XL, blue	150 pr/Cs

- *Manufactured of spunbond polypropylene*
- *Sewn seams, not glued - offers greater strength*
- *Elastic opening provides a snug comfortable fit*
- *Non-skid shoe covers offer improved floor traction*
- *Latex free*



Medicut™ Scalpels - Disposable - Sterile

Reorder No.	Description	Quantity
4110	#10	10/Bx
4111	#11	10/Bx
4112	#12	10/Bx
4115	#15	10/Bx
4120	#20	10/Bx
4121	#21	10/Bx
4122	#22	10/Bx
4123	#23	10/Bx

- *Disposable scalpels...for maximum safety, efficiency, and economy*
- *Sterile...individual peel open pouch*
- *Sheffield stainless steel blade*
- *Autoclavable plastic handle*
- *Plastic blade guard can be replaced after use*



Specimen Containers

Reorder No.	Description	Quantity
4253	4 oz sterile, tab type, tamper-evident, individually bagged	100/Cs
4254	4 oz sterile, tab type, tamper-evident	100/Cs
4256	4 oz non-sterile	Bulk 500/Cs



Cotton Tipped Applicators

Reorder No.	Description	Quantity
4301	Non-sterile 3"	10/10/100/Cs
4302	Non-sterile 6"	10/10/100/Cs
4303	Sterile 3", 2/pouch	10/100/2's/Cs
4304	Sterile 6", 1/pouch	10/100/1's/Cs
4305	Sterile 6", 2/pouch	10/100/2's/Cs

- *Cotton/Rayon tipped applicator sticks and tips are manufactured to uniform length and shape*
- *Cotton tips are highly absorbent*
- *Sterile applicators are packaged in convenient peel-down pouches*
- *Bag is autoclavable with a colored auto claving process marker, turning reddish pink to brown*

