




WE CARE LIKE FAMILY



## Endotracheal Tubes

### Used to deliver oxygen in higher concentrations

- May also be used as a route for administration of certain medications
- Help secure and maintain patient airways
- Endotracheal Tubes have a Murphy eye hole
- Sterile



### Uncuffed

Reorder No.	Description	Quantity
4530	2.0mm, Uncuffed	10/Pkg
4531	2.5mm, Uncuffed	10/Pkg
4532	3.0mm, Uncuffed	10/Pkg
4533	3.5mm, Uncuffed	10/Pkg
4534	4.0mm, Uncuffed	10/Pkg
4535	4.5mm, Uncuffed	10/Pkg
4536	5.0mm, Uncuffed	10/Pkg
4537	5.5mm, Uncuffed	10/Pkg

### Cuffed



Reorder No.	Description	Quantity
4538	4.0mm, Cuffed	10/Pkg
4539	4.5mm, Cuffed	10/Pkg
4540	5.0mm, Cuffed	10/Pkg
4541	5.5mm, Cuffed	10/Pkg
4542	6.0mm, Cuffed	10/Pkg
4543	6.5mm, Cuffed	10/Pkg
4544	7.0mm, Cuffed	10/Pkg
4545	7.5mm, Cuffed	10/Pkg
4546	8.0mm, Cuffed	10/Pkg
4547	8.5mm, Cuffed	10/Pkg
4548	9.0mm, Cuffed	10/Pkg




Endotracheal Tube with stylette on reverse side





WE CARE LIKE FAMILY

*“An intubating stylet is a malleable metal wire designed to be inserted into the endotracheal tube to make the tube conform better to the upper airway anatomy of the specific individual.”*



### Uncuffed with Stylette

Reorder No.	Description	Quantity
4549	with Stylette – 2.0mm, Uncuffed	10/Pkg
4550	with Stylette – 2.5mm, Uncuffed	10/Pkg
4551	with Stylette – 3.0mm, Uncuffed	10/Pkg
4552	with Stylette – 3.5mm, Uncuffed	10/Pkg
4553	with Stylette – 4.0mm, Uncuffed	10/Pkg
4554	with Stylette – 4.5mm, Uncuffed	10/Pkg
4555	with Stylette – 5.0mm, Uncuffed	10/Pkg



### Cuffed with Stylette

Reorder No.	Description	Quantity
4556	with Stylette – 4.0mm, Cuffed	10/Pkg
4557	with Stylette – 4.5mm, Cuffed	10/Pkg
4558	with Stylette – 5.0mm, Cuffed	10/Pkg
4559	with Stylette – 5.5mm, Cuffed	10/Pkg
4560	with Stylette – 6.0mm, Cuffed	10/Pkg
4561	with Stylette – 6.5mm, Cuffed	10/Pkg
4562	with Stylette – 7.0mm, Cuffed	10/Pkg
4563	with Stylette – 7.5mm, Cuffed	10/Pkg
4564	with Stylette – 8.0mm, Cuffed	10/Pkg
4565	with Stylette – 8.5mm, Cuffed	10/Pkg
4566	with Stylette – 9.0mm, Cuffed	10/Pkg



See our website for our complete line of Silicone Foley Catheters (No.5071-No.5089)



Find us on:



For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)