

The Brand
to Know!™




Quality Medical Products
Since 1967



Endotracheal Tube Introducer (Bougie)

Solid polypropylene stick used to assist in difficult intubations in emergency situations

- Used to guide the endotracheal tube during difficult intubations
- Tube slides over the introducer which helps to guide the way
- Especially helpful for “blind” intubations
- Sterile




Reorder No.	Description	Quantity
4582	Child, 10Fr	10/Pkg
4583	Adult, 15Fr	10/Pkg



Already mobile?
Scan your **smart phone**
here for more information.

Find us on:



See our website for our **Red Rubber Urethral & Silicone Foley Catheters**
(No.5060-No.5089)

For more information contact your
local Dynarex® distributor sales representative.
Or call 888-396-2739 to find your nearest distributor.

Dynarex • 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2015 Dynarex Corporation