



WE CARE LIKE FAMILY

## Sharps Containers

**Dynarex high quality containers are designed to be strong and lightweight**

- Containers for used medical needles and other sharp medical instruments
- Designed so needles can be dropped in the container without touching the outside
- Commonly used in EMS, hospitals, nursing homes, clinics and doctors offices, etc.
- Transparent lid allows for visual monitoring of contents to avoid overfilling

Reorder No.	Description	Quantity
4622	Sharps Container - Slide Lid, 1 qt.	60/Cs
4623	Sharps Container - Slide Lid, 2 qt.	60/Cs
4624	Sharps Container - Mailbox Lid, 5 qt.	20/Cs
4625	Sharps Container - Mailbox Lid, 5.4 qt.	20/Cs
4626	Sharps Container - Slide Lid, 1 gal.	24/Cs
4627	Sharps Container - Slide Lid, 2 gal.	24/Cs
4628	Sharps Container - Mailbox Lid, 3 gal.	12/Cs
4629	Sharps Container - Slide Lid, 6 gal.	12/Cs
4630	Sharps Shaft™, Transportable	24/Cs

\*Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.

\*Reorder No. 4624 - the 5.4qt container fits the Covidien/Kendall Wall Cabinet #85161H

\*\*Reorder No. 4625 - the 5.4qt container fits a BD locking bracket (#305447)



Also see our website for table paper (Reorder No. 4479 - 4482 & 4485 - 4487)

Already mobile? Scan your smart phone here for more information.



Find us on:



**dynarex**® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)



Tattoo

Physician



WE CARE LIKE FAMILY

Hospital



- Horizontal drop design for maximum capacity to fill line
- Each container locks for final disposal
- Nestable containers save space and reduce shipping costs
- Easy-to-see fill levels
- Single use containers which are disposed of with the waste inside



**Reorder No. 4622**

Sharps Container - Slide Lid, 1 qt. 60/Cs

Size: 4.3"×3.5"×6.5" (11×9×16.5cm)



**Reorder No. 4623**

Sharps Container - Slide Lid, 2 qt. 60/Cs

Size: 4.3"×3.5"×8.5" (11×9×21.5cm)



**Reorder No. 4624\***

Sharps Container - Mailbox Lid, 5 qt. 20/Cs

Size: 10.7"×4.6"×11.8" (27×11.5×30cm)

\*Reorder No. 4624 - the 5.4qt container fits the Covidien/Kendall Wall Cabinet #85161H



**Reorder No. 4625\*\***

Sharps Container - Mailbox Lid, 5.4 qt. 20/Cs

Size: 12"×4.5"×11.3" (30.5×11.5×28.5cm)

\*\*Reorder No. 4625 - the 5.4qt container fits a BD locking bracket (#305447)



**Reorder No. 4626**

Sharps Container - Slide Lid, 1 gal. 24/Cs

Size: 10.3"×7"×5.5" (26×18×14cm)



**Reorder No. 4627**

Sharps Container - Slide Lid, 2 gal. 24/Cs

Size: 11"×7.3"×14.2" (28×18.5×36cm)



**Reorder No. 4628**

Sharps Container - Mailbox Lid, 3 gal. 12/Cs

Size: 15"×7.3"×14.2" (38×18.5×36cm)



**Reorder No. 4629**

Sharps Container - Slide Lid, 6 gal. 12/Cs

Size: 15.2"×10"×16" (38.5×25×40cm)



**Reorder No. 4630**

Sharps Shaft™ - Transportable 24/Cs

Size: 6.37"×1.22" (16.17×3.1cm)

\*Covidien/Kendall is a registered trademark of Medtronic.®

\*\*BD is a registered trademark of Becton, Dickinson and Company.®

Visit [www.dynarex.com](http://www.dynarex.com)