



WE CARE LIKE FAMILY

Thumb Forceps



**Ideal for grasping surfaces.
Used for medical and first
aid purposes.**

Dynarex Plastic Thumb Forceps are suitable for removing splinters from the skin, and can replace the forceps found in a variety of first aid kits.

NEW!

- Tapered ends for precise pick-ups
- Sturdy and disposable
- Unique construction maintains tip alignment
- Straight serrated tip
- Non-locking
- Anti-slip thumb and finger ridges for extra control
- Not made with natural rubber latex

Reorder No.	Description	Quantity
4917	Plastic Thumb Forceps Bulk – 4.25"	3000/Cs
4918	Plastic Thumb Forceps Individually Bagged – 4.25"	1000/Cs

*Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.



See our website for our complete line of Stretch Gauze Bandages (3102-3134) and our Waterproof Adhesive Tapes in a plastic spool (3651-3657)



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com