



WE CARE LIKE FAMILY™



CPR Products

Helps Protect Rescuers In An Emergency Situation

CPR Face Shield (Reorder No. 4921)

- A lightweight face shield to prevent mouth to mouth contact during CPR
- One way valve eliminates blow back to caregiver
- Shield is individually poly bagged
- Contains detailed instructions
- Not made with natural rubber latex

CPR Face Shield (Reorder No. 4923)

- Comes in convenient soft case for easy access
- Handy swivel eye snap hook attaches to belt loop or key chain
- Not made with natural rubber latex



4921

4923

Reorder No.	Description	Quantity
4921	CPR Shield w/one way valve	100/Cs
4923	CPR Shield in soft case	100/Cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.

Visit our website for our Ammonia Inhalants (Ampule) – 33cc (Reorder No. 1401)



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2017 Dynarex Corporation

CPR Rescue Mask Kits



WE CARE LIKE FAMILY™

Makes resuscitation safer for the responder and more efficient for the patient.

Since you won't be pinching the patient's nose shut, each rescue breath will have two avenues of entry into the body. Supplemental oxygen, connected by way of the mask's oxygen inlet, increases concentrations of organ-sustaining oxygen by 50%. Remove the included valve, and the mask can be connected to a BVM. The BVM with reservoir will increase patient oxygen concentrations to 100%.

CPR Rescue Mask Kit (Reorder No. 4924)

- Completely seals to the faces of infants, children and adults
- Includes 1 pair of vinyl exam gloves
- Contains detailed instructions
- CPR Mask with optional oxygen port
- Hard case with wrist strap
- Not made with natural rubber latex



Reorder No.	Description	Quantity
4924	CPR Rescue Mask Kit	25/Cs

Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.

Visit our website for our **Ammonia Inhalants (Ampule) – 33cc** (Reorder No. 1401)



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2017 Dynarex Corporation