



## Urethral Catheterization Tray



### Trust in Dynarex for Patient Care

- **Straight tip**
- **Latex-free**
- **Medical grade PVC**
- **Offset polished eyelets to maximize urine flow**



#### Contents:

- (1) 14 French (male) intermittent catheter
- (1) 1000cc outer basin tray
- (1) Pair of vinyl gloves
- (1) Waterproof drape (blue, 45 x 42 cm)
- (1) Fenestrated drape (white, 45 x 42 cm)
- (1) Lubricating jelly packet (3g)
- (3) Povidone Iodine swabsticks
- (1) 4 oz specimen cup
- (1) Patient label



Reorder No.	Urethral Catheterization Tray	Quantity
4928	Latex Free, Sterile	20/cs

See our website for our complete range of intermittent catheters (#4981-4986, #4991-4995)

