

The Brand  
to Know!™



Quality Medical Products  
Since 1967



5605

## Nebulizers

**For treatment of asthma,  
allergies and other  
respiratory disorders**

**Converts prescribed medication into  
aerosol mist for easy inhalation.**

Reorder No.	Description	Quantity
5605	Elite Compressor Nebulizer	4/Cs
5606	Portable Compressor Nebulizer	6/Cs

5605-5606 HCPCS: E0570

**Now Comes With a  
Three Year Warranty!**

### Elite Compressor Nebulizer 5605

- Includes standard accessories: air tube, mouthpiece, filters (5 pcs), adult mask, child mask.
- All accessories can be stored in convenient compartment.
- Power: 120V-60Hz (Standard U.S. 110/120V – 50 or 60 Hz).
- Warrantied for three years – power cord 2 years. Must complete warranty registration card.
- Compare to Drive Medical AeroNeb®\* (#7090-7091).
- Compressor Pressure Range is 30-50 PSI.

### Portable Compressor Nebulizer 5606

- Ultra compact design with handle, making it extremely portable
- Includes standard accessories: air tube, mouthpiece, filters (5 pcs), adult mask, child mask.
- Power: 120V-60Hz (Standard U.S. 110/120V – 50 or 60 Hz).
- Warrantied for three years – power cord 2 years. Must complete warranty registration card.
- Compare to Probasics Rite-Neb 3™\*\* – But includes mask, which they don't.
- Compressor Pressure Range is 30-50 PSI.

\*Drive Medical AeroNeb is a registered trademark of Aerogen, Inc.

\*\*Probasics Rite-Neb 3 and the design are trademarks of Garden City Medical, Inc.



5606

Nebulizer Standard Filter Replacement Information for item #5605-5606		
Item #	Description	Standard Replacement Filters - Available on-line reference #:
5605	Elite Compressor Nebulizer	Item #: PMI805SLTI
5606	Portable Compressor Nebulizer	Item #: PMI805SLTI

For more information contact your  
local Dynarex® distributor sales representative.  
Or call 888-396-2739 to find your nearest distributor.