

The Brand
to Know!™



Quality Medical Products
Since 1967

Sterile Nitrile Exam Gloves



Exam-Blue
Pairs

Exam-Blue
Singles

**Ideal for use in OB/GYN,
ER/Trauma, Burn Centers,
and Pharmacy, or other
areas where durable latex-
free protection is needed.**

- Powder free for easier donning
- Ambidextrous
- Textured fingertips for improved grip
- Beaded cuff protects wrists and arms from messy runoff (Cuff size is 9")
- Sterile
- Packaged in inner wallet (single or pairs)
- Glove thickness is 4 mil.

See our website for our
Sterile Nitrile Surgical Gloves
(6533-6539)

Reorder No.	Description	Quantity
6522	Exam Gloves-S, Singles-Blue	8/100/Cs
6523	Exam Gloves-M, Singles-Blue	8/100/Cs
6524	Exam Gloves-L, Singles-Blue	8/100/Cs
6528	Exam Gloves-S, Pairs-Blue	8/50/Cs
6529	Exam Gloves-M, Pairs-Blue	8/50/Cs
6530	Exam Gloves-L, Pairs-Blue	8/50/Cs



Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.

Find us on:



Scan your
smart phone
here for more
information.



For more information contact your
local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

Dynarex • 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2015 Dynarex Corporation