



WE CARE LIKE FAMILY™



EpiRite™ Syringe



Easy to Use and Convenient

- 1 mL Luer Slip
- Easy to read markings - 0.1 mL, 0.15 mL, and 0.3 mL
- Allows for fast and easy dosing
- Syringe provided empty, single use only



Reorder No.	Product Name	Quantity
6985	EpiRite Syringe (1 mL Luer Slip)	10/100/Cs



Find us on:



dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2019 Dynarex Corporation

R190214