



WE CARE LIKE FAMILY



IV Administration Sets and IV Extension Sets

Available in various lengths and dimensions to accommodate various applications

- Non-DEHP tubing
- Not made with natural rubber latex
- Non-Pyrogenic
- All IV sets are gravity sets.
Not compatible with infusion pump.
- Sterile-individually packaged
- Sold by inner box



IV Administration Sets

Reorder No.	Description	Quantity
7030	10 Drop/mL, 83"	50/Box
7031	15 Drop/mL, 83"	50/Box
7032	15 Drop/mL, 78"	50/Box
7033	15 Drop/mL, 106"	50/Box
7034	15 Drop/mL, 73"	50/Box
7035	15 Drop/mL, 84"	50/Box
7036	15 Drop/mL, 85"	50/Box
7037	15 Drop/mL, 92"	50/Box
7048	20 Drop/mL, 100"	50/Box
7049	20 Drop/mL, 92"	50/Box
7050	20 Drop/mL, 92"	50/Box
7051	60 Drop/mL, 83"	50/Box
7052	60 Drop/mL, 78"	50/Box
7053	60 Drop/mL, 75"	50/Box

IV Extension Sets


Reorder No.	Description	Quantity
7060	6", Needle-Free Luer Lock Connector	100/Box
7061	7", Needle-Free Luer Lock Connector	100/Box
7062	8", Needle-Free Luer Lock Connector	100/Box

Compares to Amsino, B. Braun, Exel and MedSource*.

*© Amsino International. © B. Braun Medical Inc. © EXEL International. © MedSource.



See our website for our full line of syringes and needles (#6900-7015)



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.


10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2016 Dynarex Corporation

dynarex IV Administration Sets and IV Extension Sets



See our website for
all our medical products!
www.dynarex.com




IV Administration Sets

Reorder No.	Drip-Rate and Length	Quantity	Ports & Port Type (all rotating Luer Lock connector)	Quantity of Luer Lock Connectors (needleless)	Quantity of Y-Injector Sites (needle)	Quantity of Micron Filters	Quantity of Accurate Regulators
7030	10 Drop/mL, 83"	50/Box	1 Y-injection site and 1 Luer Lock connector (needleless*)	1	1	-	-
7031	15 Drop/mL, 83"	50/Box	1 Y-injection site and 1 Luer Lock connector (needleless*)	1	1	-	-
7032	15 Drop/mL, 78"	50/Box	1 Y-injection site	-	1	-	-
7033	15 Drop/mL, 106"	50/Box	2 Y-injection sites	-	2	-	-
7034	15 Drop/mL, 73"	50/Box	1 Y-injection site	-	1	-	-
7035	15 Drop/mL, 84"	50/Box	2 Luer Lock connectors (needleless*)	2	-	-	-
7036	15 Drop/mL, 85"	50/Box	Without injection site	-	-	-	-
7037	15 Drop/mL, 92"	50/Box	2 Y-injection sites	-	2	-	-
7048	20 Drop/mL, 100"	50/Box	1 Y-injection site w/micron filter	-	1	1	-
7049	20 Drop/mL, 92"	50/Box	1 Y-injection site w/micron filter and regulator	-	1	1	1
7050	20 Drop/mL, 92"	50/Box	1 Y-injection site w/0.2 micron filter	-	1	1	-
7051	60 Drop/mL, 83"	50/Box	1 Y-injection site and 1 Luer Lock connector (needleless*)	1	1	-	-
7052	60 Drop/mL, 78"	50/Box	2 Y-injection sites	-	2	-	-
7053	60 Drop/mL, 75"	50/Box	Without injection site	-	-	-	-

IV Extension Sets

Reorder No.	Description	Quantity
7060	6", Needle-Free* Luer Lock Connector	100/Box
7061	7", Needle-Free* Luer Lock Connector	100/Box
7062	8", Needle-Free* Luer Lock Connector	100/Box



Needle-free Luer Lock connector for injection.

*All needleless and needle-free products include the Needleless PRN Rotary Luer Lock Connector



Needleless PRN Connector also available separately, Reorder No. 7054.

Coming Soon!

dynarex® is your one-stop medical supply company since 1967



Find us on:

