



WE CARE LIKE FAMILY™



IV Administration Sets and IV Extension Sets

Available in various lengths and dimensions to accommodate various applications

- Non-DEHP tubing
- Not made with natural rubber latex
- Non-Pyrogenic
- All IV sets are gravity sets. Not compatible with infusion pump
- Sold by the inner box



IV Administration Sets

Reorder No.	Description	Quantity
7030	10 Drop/mL, 83"	50/Bx
7031	15 Drop/mL, 83"	50/Bx
7032	15 Drop/mL, 83"	50/Bx
7033	15 Drop/mL, 83"	50/Bx
7034	15 Drop/mL, 83"	50/Bx
7035	15 Drop/mL, 83"	50/Bx
7036	15 Drop/mL, 83"	50/Bx
7037	15 Drop/mL, 83"	50/Bx
7038	10 Drop/mL, 37" w/ Hanger & Lever Lock Cannula	50/Bx
7039	10 Drop/mL, 83", 1 Injection Site	50/Bx
7040	10 Drop/mL, 83" w/ 6" extension, 1 Injection site, 1 Needleless Port	50/Bx
7041	10 Drop/mL, 89", 2 Injection Sites, Check Valve	50/Bx
7042	10 Drop/mL, 112", 2 Needleless Ports, Check Valve	50/Bx
7043	15 Drop/mL, 106", 2 Needleless Ports	50/Bx
7044	20 Drop/mL, 92", 1 Injection Site	50/Bx
7045	20 Drop/mL w/ Regulator, 92", 1 Needleless Port w/ .02 Micron Filter	50/Bx
7046	20 Drop/mL w/ Regulator, 100", 1 Needleless Port w/ 1.2 Micron Filter	50/Bx
7047	20 Drop/mL, 100", 1 Needleless Port w/ .02 Micron Filter	50/Bx
7048	20 Drop/mL, 100"	50/Bx
7049	20 Drop/mL, 92"	50/Bx
7050	20 Drop/mL, 92"	50/Bx
7051	60 Drop/mL, 83"	50/Bx
7052	60 Drop/mL, 78"	50/Bx
7053	60 Drop/mL, 75"	50/Bx

NEW Products!

IV Extension Sets

Reorder No.	Description	Quantity
7060	6", Needle-Free Luer Lock Connector	100/Bx
7061	7", Needle-Free Luer Lock Connector	100/Bx
7062	8", Needle-Free Luer Lock Connector	100/Bx
7064	IV Extension Set - 18" w/ Regulator	50/Bx
7065	IV Extension Set - 36"	50/Bx

dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2018 Dynarex Corporation

R180803



Find us on:

