



WE CARE LIKE FAMILY

Finger Pulse Oximeters

Non-Invasive Oxygen Saturation Measuring

Finger Pulse Oximeter

Calculates how much oxygen is carried in your blood and monitors saturation of a patient's hemoglobin. It is also a convenient non-invasive measurement method to detect hypoxia.

- Unbeatable accuracy – clinically tested
- Drop tested durable and moisture proof
- Perfect spot checks or continuous readings
- Simple one-button operation
- Checks heart rate
- Monitors respiratory function
- Monitors a patient's oxygenation and pulse rate when undergoing surgical procedure
- Comes individually in a carton with instructions
- Requires two AAA batteries (not included)
- One-year limited warranty – warranty card must be completed
- Suitable for finger thickness .37" - .86" (9.5mm - 21.9mm)
- Not made with natural rubber latex



7088

NEW!



7090

Reorder No.	Description	Quantity
7088	Finger Pulse Oximeter	10/Cs
7090	Pediatric Finger Pulse Oximeter	10/Cs

NEW!

Pediatric Finger Pulse Oximeter

A portable non-invasive device intended for spot checking of oxygen saturation of arterial hemoglobin and pulse rate of patients at home and in hospital.

- Child-friendly instrument face
- Specially marked to tell children where to place their fingers
- Six changeable display modes
- Reads SpO₂ and pulse rate
- Requires two AAA batteries (not included)
- One-year limited warranty – warranty card must be completed
- Suitable for finger thickness .2" - .75" (5mm - 19mm)
- Not made with natural rubber latex



See our website for other diagnostic equipment

Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.

dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com



Find us on:

