Aluminum Mechanical Height Rod

KEY FEATURES

Precision engineered and constructed of heavy gauge aluminum, the Doran DS1150 Mechanical Height Rod will deliver accurate height readings for years and years. Flip the measurement arm up and slide the height measurement arm into position. The patient's height is highlighted with an easy-to-read, bright red line.

Precisely marked in centimeters and inches, the DS1150 can be used to measure adult or pediatric patients.

The DS1150 is also very easy to install on our DS7060 and DS7100 handrail scales. When not in use, the rod lays flat – so it's ideal for mounting against the wall or conveniently behind a door.

SPECIFICATIONS

Warranty:	Two Years
Measuring Range:	24 - 82" (60 - 210 cm)
Graduations:	1/16" (1 mm)
Construction:	Aluminum and ABS Plastic
Dimensions:	Height Rod 59" L x 3" W x 1" H
	Measurement Arm 11" L x 8" W x 3.5" H





Doran's height rod features our signature red markings that ensure accurate and easy-to-read height measurements.

35