



EarthLite®
WORLD'S #1 BRAND IN MASSAGE™

LUNA™
Simple yet Elegant Maple Treatment Tables



HEALING HUMANITY THROUGH TOUCH™

LUNA™

Full Sized, Lightweight Massage Table

The Luna is an innovative Reiki-ready aluminum frame and our Pro-Lite™ cushioning system combine to create a super-lite, yet comfortable portable table. Add the Shimmy Lok™ leg brace system and you have arguably the lightest, strongest, most stable and comfortable full sized massage table on the market today. Available as table only or in a value package.

Reiki-Ready Super-Lite Frame

- Expertly designed by Earthlite's advanced engineering team
- Only 29 pounds – easy to transport
- Greater client access for specific massage work
- Innovative dual post leg stabilizers

Innovative Components

- Patent pending Shimmy-Lok™ Leg Brace System eliminates leg instability inherent in other aluminum frame designs
- Patented Cradle-Lock™ cabling system with 1000 lb. rated aircraft grade steel cables

Luxurious Upholstery

- Pro-Lite™ Cushioning System
- Buttery soft, yet eco-friendly Natursoft™ upholstery

Available as Value Package

Includes Flex-Rest™ Self Adjusting Facecradle with FacePillow and durable carry case.



Features

- Patent pending Reiki-ready super-lite aluminum frame
- Full-length piano hinge
- Patent pending Shimmy-Lok™ leg brace system
- Mid-brace™ cradle-lock cabling system with 1,000 lbs test cables
- Pro-Lite™ 2-layer 2-¼ cushioning system
- Natursoft™ Upholstery – 100% PU
- Dual Maple headrest outlets

Specifications

- Width: 30" (76 cm)
- Length: 73" (185 cm)
- Height Range: 24" - 34" (61 cm - 86 cm)
- Static Weight: 3,000 lbs (1,361 kg)
- Max. Working Weight: 750 lbs (340 kg)
- Lifetime Limited Warranty
- CE



Flex-Rest™ Self-Adjusting Platform

Flexes and pivots to adapt to the unique shape of your client's face.



Natursoft™ Upholstery

