



EarthLite®  
WORLD'S #1 BRAND IN MASSAGE™

## ELLORA™

THE BEST VALUE IN SPA

The Ellora™ line is by far the most popular lift table in the market, it offers a lifetime of use at an affordable price. This heavy duty ADA Compliant table features a solid steel frame with a reliable motor that holds up to 600 lbs. Earthlite's Ellora™ provides you with the perfect balance between reliability and comfort, with deluxe cushioning and upholstery it provides years of trouble-free use.



HEALING HUMANITY THROUGH TOUCH™

# ELLORA COLLECTION

## Features

- Heavy duty steel frame
- Reliable actuator
- Pro-Plush™ Deluxe 3", 3-layer cushioning system
- 100% PU Natursoft™ upholstery in 6 standard colors: Black, Mystic Blue, Latte, Sterling, Marie's Beige, Vanilla Crème
- Hands free foot control
- Wheels for ease of movement
- ADA compliant
- Lifetime Limited Warranty
- 3 years Warranty on foam and upholstery
- 2 years Warranty on electric parts

## Options

- Manual tilt top
- Non-electric pneumatic tilt or salon top
- Electric tilt or salon top with hand control
- Extra foot pedal
- Additional Natursoft™ color selections
- Ultraleather® or Promessa® upholstery
- Flex arms, neck roll and salon pillow
- Flex-Rest™ or Caress™ platform
- Conforma LE™ Breast Recess Comfort System
- 4" Memory foam system
- ETL option

## Specifications

- Widths: 28", 30", 32"
- Length: 73"
- Height range: 18" – 37"
- Recommended lift capacity: 600 lbs.
- CE



Flat




Tilt



Salon


### Pro-Plush™ Deluxe 3-Layer, 3" Cushioning



### Non-Electric Pneumatic Adjustment



### Manual Tilt Adjustment



### Conforma LE™ Breast Recess Comfort System



The system eliminates breast discomfort and improves spinal alignment for women and men during face down spa and medical treatments.

Clients can make discrete adjustments with easy to reach and intuitive controls.

