



EarthLite®  
WORLD'S #1 BRAND IN MASSAGE™

APEX™  
STATIONARY TABLE




HEALING HUMANITY THROUGH TOUCH™

# APEX™ STATIONARY TABLE

If you must own the very best, our Apex Stationary Bench is for you. It is the ultimate adjusting bench with tilting headpiece. Our Apex features 3 inches of specialized foam, selected specifically for patient comfort and doctor functionality. You'll appreciate the Apex heavy-duty metal construction and superior foam system with extra cushioning on the headpiece. Its modern design and seamless upholstery makes our Apex a true "Doctor's Choice".

## Dimensions


- Length: 72" (183 cm)
- Width: 24" (61 cm)
- Height: 19" to 28" (48 cm - 71 cm)
- Weight: 101 pounds (46 kg)
- Warranty: Lifetime on frame, 3 years on vinyl covering, 3 years on foam



## Features

- Modern design
- Seamless upholstery
- 3 inch high-density foam
- Paper roll holder and cutter attachment
- Stationary armrest
- Powder coated white steel frame

## Nature's Touch Upholstery Options



Tilting headpiece with 4" of cushioning.



Manually adjustable height.



Powder coated steel frame.