


# EU25 TRISTAR

3 Section Therapeutic Table (No Drop Foot)



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue   Licorice Black   Dim Grey   Merlot (Red Wine)   Creamy Beige



1 Motor



Hand Control




Foot Control



2 / 3 / 4 Sides Elevating Bar

**Maximum User Weight**  
225 kgs / 496 lbs



## Product Features

This is an EU25 - 3 Section Therapeutic Table with No Drop Foot. The leg section is further supported by a bracket frame underneath to make this table a 3 Section with "No Drop Foot." The EU25 has the following features:

- Hydraulic version (EU25H) is also available
- Short Head Section (Length: 41cm / 16.1" inch)
- Long Leg Section (Length: 85cm / 33.5" inch)
- Equipped with 1 Motorized Actuator for Hi/Lo Adj.
- Twin "Oversized" Pillars for Added Stability
- Individual "Leveling Foot" on All 4 Corners
- Retractable Castors at Both Ends
- Twin Gas Struts for Solid Head and Leg Sections Support
- Bracket Frame with 2 Welding Points for Solid Support

## Technical Specifications

