

EU300

TRAÇION-t

4 Section Traction Table (Electric)

EU300 Traction is the ultimate Traction Table for all Traction Needs: Lumbar + Cervical Traction



- 1 Year on Gas Struts, Castors, Upholstery
- 2 Years on Electronic Components
- 3 Years on Frame & Pillar



Prussian Blue Licorice Black Dim Grey Merlot (Red Wine) Creamy Beige



2 Motors



Hand Control



Foot Control



2 / 3 Sides Elevating Bar

Maximum User Weight
225 kgs / 496 lbs

Product Features

ALL FOUR CORNERS OF THE EU300 ARE SUPPORTED BY 3 INCH SWIVEL CASTORS, ONE AT EACH CORNER.

- 1st Motorized Actuator for Hi/Lo Adj.
- 2nd Motorized Actuator for Centre Arching Mechanism
- Twin “Oversized” Pillars for Added Stability
- Individual “Leveling Foot” on All 4 Corners
- Retractable Castors at Both Ends
- Twin Gas Struts for Solid Head Section Support
- Maximum Load Capacity of 225 kgs / 496 lbs



Special Features

The EU300 Taction is equipped with ALL NEW “Centre Arching Mechanism” for Motorized LUMBAR ELEVATION. The facilitates any pulling traction force on the lumbar area with a slight upward lifting support (up to 26° RISE) especially when the patient is lying in the supine position.



Dual TSS Switches

For easy accessibility & control, the “Centre Arching/Lifting” is controlled by 2x TSS Switches on both sides just underneath the mid-section of table top.



Adjustable Recoil Spring

The “Adjustable Recoil Spring” helps the 3rd and 4th section of table top gradually returns to its original position as the applied traction force decreases. We implemented a handle at the end of of the table (near the 4th section) for recoil spring intensity adjustment. With this feature, therapists can adjust the intensity anytime during the therapy to make the decompression smooth and effective, depending on the patient’s weight.