



DNV BUSINESS ASSURANCE

EC CERTIFICATE - FULL QUALITY ASSURANCE SYSTEM

Certificate No. 92630-2011-CE-USA-NA 2.0 This Certificate consists of 3 pages

This is to certify that the Quality Management System of

EXELINT International Co.

2500 Santa Fe Avenue, Redondo Beach, CA 90278, USA for design, production and final product inspection/testing of

Disposable Medical Devices

has been assessed with respect to the conformity assessment procedure described in Article 11.1 a and annex II and Article 11.3 a and Annex II excluding section 4 (Module H) of Council Directive 93/42/EEC on Medical Devices, as amended, and found to comply

Further details are given overleaf ain.

Place and date: Høvik, 13 November 2014

Aud Løken Eiklid

30 January 2019

For DET NORSKE VERITAS CERTIFICATION AS NORWAY

NORWEGIAN ACCREDITATION PROD 002

Jenny Helen Nytun Technical Reviewer

EC Design-Examination Certificate Directive 93/42/EEC on Medical Devices, Annex II, section 4



Design Approval no. 92630-2011-CE-USA-NA-D

Manufacturer address:		
500 Santa Fe Avenue		
Redondo Beach, CA 90278 - USA		
ype of medical device and identification no.:	Class of Medical Device:	GMDN code:
Anaesthesia (Epidural, Spinal) Needle	ш	12729
thort description of the medical device:		
The device is EO sterilized and used to inject anaesthetic	es into spine or extract spinal fl	uid from the
pine.		

Limitations:

Any changes in the Design shall immediately be reported to Det Norske Veritas Certification AS in order to examine whether this Certificate remains valid. Annual Periodical Audits will be held to verify the validity of this Certificate.

his certificate is valid until: 30 January 2019 for Det Norske Veritas Certification AS



Net -

CERTIFICATE



This is to cartify that the company

EXELINT International, Co.

ISO 13485 : 2003



4. Shelit



Comfort point INSULIN SYRINGES

0.3cc LO-DOSE INSULIN SYRINGE

with permanently attached needle

Catalog No. Description

26018 U-100 29G x 1/2" 26014 U-100 30G x 5/16" 26010 U-100 31G x 5/16"

10/bag,100/box,500/case

1/2cc LO-DOSE INSULIN SYRINGE

with permanently attached needle

Catalog No. Description

26026 U-100 28G x 1/2" 26028 U-100 29G x 1/2" 26015 U-100 30G x 5/16" 26011 U-100 31G x 5/16" 10/bag,100/box,500/case

1cc INSULIN SYRINGE

with permanently attached needle

Catalog No. Description

26027 U-100 28G x 1/2" 26029 U-100 29G x 1/2" 26016 U-100 30G x 5/16" 26012 U-100 31G x 5/16"

10/bag,100/box,500/case

 This product does not contain natural rubber latex.

- Sterile
- Non-toxic
- Non pyrogenic
- Permanently attached needle
- Reduced dead space for reduced Insulin waste
- Smooth plunger action for greater filling accuracy
- Auto ground needle tip for greater comfort





Comfort point INSULIN PEN NEEDLES


- Optimum treatment with maximum simplicity
- This device does not contain natural rubber Latex.
- Sterile
- Non-toxic
- Non-pyrogenic
- Double point system
- Patient point with triple beveling for easy and atraumatic penetration
- Cartridge point sharpened for a perfect punch with no fragmentation
- Compatible with all available pen models

Easy flow of medication without any

INSULIN PEN NEEDLES

Catalog No.	Description
26003	29 G (0.33) x 1/2"(12mm)
26004	31 G (0.25) x 1/4"(6mm)
26005	31 G (0.25) x 1/3"(8mm)
26006	31 G (0.25) x 1/6"(4mm)

100pcs/box, 12box/case (1200pcs/case)



TUBERCULIN SYRINGES ALLERGY SYRINGES

1cc TUBERCULIN SYRINGE

with detachable needle, Luer-slip

Catalog No. Description

26040* 1cc 27Gx1/2" 26042* 1cc 25Gx1" 26043* 1cc 26Gx1/2" 26044* 1cc 25Gx5/8"

100/box,1,000/case

1cc TUBERCULIN SYRINGE ONLY

Catalog No. Description

26048* 1cc Slip Tip with cap 26049* 1cc Lock Tip with cap

26050 1cc Lock Tip 100/box,1,000/case

1cc ZERO DEAD SPACE TUBERCULIN SYRINGE

with permanently attached needle

Catalog No. Description

26046 1cc 25G x 5/8" 10/bag,100/box, 500/case

1cc ALLERGY SYRINGES

Catalog No. Description

26061 1cc 27G x 1/2" 26063 1cc 27G x 3/8"

25/tray, 40/trays/case



- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Non pyrogenic
- Reduced dead space for reduced medication waste
- Smooth plunger action for greater filling accuracy
- Auto ground needle tip for greater comfort



3CC SYRINGES







* 3cc SYRINGE only

Catalog No. Description

26200 Luer Lock with cap 26201 Luer Slip with cap 100/box, 1,000/case

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Pyrogen free
- Low Dead Space
- With Cap
- Positive plunger stop
- Crystal clear barrels
- Easy to read large, bold print

* LUER -LOCK TIP

3cc Syringe/Needle Combination

Catalog No.		lo. Description	Hub Color
	26100	25G x 5/8"	Orange
	26101	23G x 1"	Blue
	26102	22G x 1"	■ Black
	26103	22G x 1 1/4"	■ Black
	26104	22G x 1 1/2"	■ Black
	26105	21G x 1 "	Green
	26107	21G x 1 1/2"	Green
	26108	20G x 1 "	Yellow
	26109	20G x 1 1/2"	Yellow
	26110	18G x 1 1/2"	Pink
	26111	25G x 1 "	Orange
	26112	25G x 1 1/2"	Orange
	26115	22G x 3/4"	■ Black
	26117	23G x 1 1/2"	Blue
	26127	21G x 1 1/4"	Green
	26129	27G x 1 1/4"	Medium Grey
	100/box,	1,000/case	

* LUER -SLIP TIP

2

3

ml

3cc Syringe/Needle Combination

Catalog N	o. Description	Hub Color
26601	22G x 3/4"	■ Black
26602	25G x 5/8"	Orange
26603	23G x 1"	Blue
26605	22G x 1 "	■ Black
26608	21G x 1 "	Green
26611	20G x 1 "	Yellow
100/box,	1,000/case	

5-6CC SYRINGES 10-12CC SYRINGES CURVED TIP SYRINGES

5-6cc Syringe/ Needle Combination LUER -LOCK TIP

Catalog No.	Description	Hub Colo
26210	22G x 1"	Black
26211	22G x 1 ½"	Black
26212	21G x 1"	Green
26213	21G x 1 ½"	Green
26214	20G x 1"	Yellow
26218	21G x 1 1/4"	Green
26219	22G x 1 1/4 "	Black
100/box, 800	/case	

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Pyrogen free
- Positive plunger stop to prevent spillage when drawing
- Smooth plunger motion

5-6cc SYRINGE only

Catalog No. Description

26230 Luer Lock with cap 26231 Luer Slip with cap

100/box, 800/case



10-12cc Syringe/ Needle Combination LUER -LOCK TIP

Catalog No.	Description	Hub Color
26250	22G x 1"	Black
26251	22G x 1 ½"	Black
26252	21G x 1"	Green
26253	21G x 1 ½"	Green
26254	20G x 1"	Yellow
26255	20G x 1 ½"	Yellow
26257	18G x 1 "	Pink
100/box, 800	/case	

10-12cc SYRINGE only

Catalog No. Description

26265 Luer Lock with cap 26266 Luer Slip with cap

100/box, 800/case

CURVED TIP SYRINGE

Catalog No. Description

26267 Curved Tip

50/box, 500/case

20-25CC SYRINGES 30-35CC SYRINGES 50-60CC SYRINGES CULINARY SYRINGES

without needle

20-25cc SYRINGE only

Catalog No. Description

26280 Luer Lock with cap 26281 Luer Slip Ecc w/cap

50/box, 500/case

30-35cc SYRINGE only

Catalog No. Description

26290 Luer Lock with cap 26291 Luer Slip Ecc with cap

26292 Cath Tip Ecc

50/box, 250/case

50-60 SYRINGE only

Catalog No. Description

26300 Luer Lock with cap 26301 Luer Slip Ecc with cap 26304 Catheter Tip with cap

25/box, 150/case

CULINARY SYRINGE

Catalog No.Description

26140 Culinary Syringe

with Culinary needle

100/bag



- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Pyrogen free
- Positive plunger stop to prevent spillage when drawing
- Smooth plunger motion





BULB SYRINGES

Catalog No. Description 25531 3 oz



HYPODERMIC NEEDLES

Regular Bevel

Catalog No.	Description	Hub Color
26400	27G x 1/2"	Medium Grey
26402	26G x 1/2"	Brown
26403	25G x 5/8"	Orange
26405	25G x 1"	Orange
26406	25G x 1 1/2"	Orange
26407	23G x 3/4"	Blue
26408	23G x 1"	Blue
26409	23G x 1 1/2"	Blue
26410	22G x 3/4"	■ Black
26411	22G x 1"	■ Black
26412	22G x 1 1/2"	■ Black
26414	21G x 1"	Green
26416	21G x 1 ½"	Green
26417	20G x 1"	Yellow
26418	20G x 1 1/2"	Yellow
26419	18G x 1"	Pink
26420	18G x 1 1/2"	Pink
26424	19G x 1"	Cream
26426	27G x 1 1/2"	Medium Grey
26427	27G x 1 1/4"	Medium Grey
26430	20G x 3/4"	Yellow
26431	20G x 1/2"	Yellow
26438	19G x 1 1/2"	Cream

100/box, 2,000/case

ALUMINUM HUB BLUNT NEEDLES

Catalog No.	Description
26471	15G x 1 1/2"
26472	17G x 1 1/2"
26473	18G x 1"
26474	19G x 1 1/2"
26475	20G x 1 1/2"
25/box,100/c	ase

NEEDLE LENGTH								
	Inches	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2 "
	Millimeters	13mm	16mm	19mm	25mm	32mm	38mm	50mm




SPECIALITY USE HYPODERMIC NEEDLES

Catalog No.	Description	Hub Color
26437	30G x 1/2"	Yellow
26439	30G x 1"	Yellow
100/box, 2,0	000/case	
26448	28G x 3/4"	Blue Green
100/box, 1,0	000/case	
26391	31G x 1/2"	Orange
26392	32G x ½"	Green
26392S	32G x ½"	Green Soft Pack
26393	33G x 1/2"	☐ Clear
100/box/10	00/case	

- This product does not contain natural rubber latex.
- Sterile, Non-toxic, Pyrogen free
- SUS 304 stainless steel
- Needles and boxes are color-coded for easy identification
- Each needle is individually packaged, sterilized and ready for use
 - Sterile, single-use, individually packaged blunt needle
 - Stainless steel laser welded, non-lubricated cannula
 - Aluminium hub, luer lock design
 - Hub crimp cannula attachment
 - Color-coded for easy identification
 - 2-piece polypropylene needle cap and sheath

SAFETY HYPODERMIC NEEDLES






Secure Iouch

SAFETY HYPODERMIC NEEDLES


Catalog No.	Description
27403	25G X 5/8"
27405	25G X 1"
27406	25G X 1 1/2"
27417	20G X 1"
27420	18G X 1 1/2"
27426	27G X 1 1/2"
27427	27G X 1 1/4"
27437	30G X 1/2"

100/box, 1000/case

This device does not contain natural rubber latex.



OR



Secure louch

SAFETY NEEDLES

- One-Handed Activation
- No solid surface required.
- Simple finger stroke activation
- Protective arm keeps fingers behind the needle.
- An audible "click" as well as tactile feel and a visual indicator confirm safety mechanism engagement.
- Needle remains locked inside activated cover.
- Fits any brand syringe.
- Safety shield can be pushed back comparative to the needle, allowing a clear view of the injection site.



As shown on a 3CC Luer Lock Syringe

This device does not contain natural rubber latex.

BLOOD COLLECTION SETS SAFETY BLOOD COLLECTION SETS



SAFETY BLOOD COLLECTION SET

Catalog No. Description

Tube **Hub Color** SAFETY PSV/M-S,L-A 21G x 3/4" 12" Green SAFETY PSV/M-S.L-A 23G x 3/4" 12" Light Blue SAFETY PSV/M-S,L-A 25G x 3/4"

12"

Dark Blue

50/box 200/case

27764

27766

27768

This device does not contain







This device does not contain natural rubber latex.

Catalog No.	Description	Tube	Hub Color		
26764	21G x 3/4"	12"	Deep Green		
26766	23G x 3/4"	12"	Deep Blue		
26768	25G x 3/4"	12"	Orange		
50/box 200/case					

SAFETY SCALP VEIN SETS SCALP VEIN SETS



Securelouch

- This device does not contain natural rubber Latex
- Individually packaged
- Thin wall stainless steel cannula, siliconized for easy insertion
- Flexible wings hold the needle in position and help to prevent inadvertent pull out
- Color-coded for easier identification



SAFETY SCALP VEIN SET

SAFETY PSV / BUTTERFLY / SMALL VEIN INFUSION

Catalog N	10.	Descr	iption
	Tube	Hub C	Color
27704	21G x 3/4"	12"	Green
27706	23G x 3/4"	12"	Light Blue
27708	25G x 3/4"	12"	Blue
50/box,50	00/case		



SCALP VEIN SET

PSV / BUTTERFLY / SMALL VEIN INFUSION

Catlog	No. Description	Tube	Hub Color
26702	19G x 3/4"	12"	Cream
26703	20G x 3/4"	12"	Yellow
26704	21G x 3/4"	12"	Deep Green
26706	23G x 3/4"	12"	Deep Blue
26708	25G x 3/4"	12"	Orange
26709	27G x 3/4"	12"	Medium Grev

50/box,500/case

AV FISTULAS - SINGLE PACK

- This product does not contain natural rubber latex.
- Ultra-thin wall needles allow maximum blood flow
- Ultra-sharp for easier penetration
- Siliconized for smoother, more comfortable penetration
- Color coded wings for easier identification
- Easy removal of luer cap



AV FISTULA - SINGLE PACK with BACKEYE
12" TUBE, CLAMP, FIXED WING

Catalog No. Description

26786 16G x 1" 26791 17G x 1"

50pcs/box,500pcs/case



IV ADMINISTRATION SETS FILTERED IV ADMINISTRATION SETS



I.V. ADMINISTRATION SET VENTED/NON VENTED

With Y-injection Site, Vented / Non Vented, With Clamp ● 78" Tubing

Priming Volume 23ml/Approx Inside Diameter 3.0 ± 0.05mm Outside Diameter 4.1 + 0.10 / -0.05



Catalog No.Description

29081 15 Drop/ml Luer Slip29083 60 Drop/ml Luer Slip

50 pcs/case



- This product does not contain natural rubber latex.
- Sterile
- Non toxic
- Non pyrogenic fluid path
- Vented/ Non Vented
- Roller clamp
- Occlusion clamp
- Pouch packaging

Catalog No.Description

29082 15 Drop/ml Luer Lock 29085 60 Drop/ml Luer Lock

50 pcs/case



Catalog No.Description

29086 15 Drop/ml Lock/Slip

27086 15 Drop/ml Lock/Slip with Needle-Free Injection Site

50 pcs/case



Catalog No.Description

29088 15 Drop/ml Lock/Slip, Vented/Non Vented, 2 Injection Sites,

Roller Clamp, 2 Occlusion Clamps ● 105" Priming Volume 28ml/Approx

50 pcs/case

DEHP free tubing

- 0.2 Micron IV filter
- Filters are constructed with materials that pass USP Class VI

FILTERED I.V. ADMINISTRATION SET VENTED/NON VENTED



Catalog No.Description

26085 20 Drop/ml Luer Lock, Vented/Non Vented, Latex-Free Injection Site,

Roller Clamp, Occlusion Clamp ● 84" Tubing

27085 20 Drop/ml Luer Lock, Vented/Non Vented, with Needle-Free Injection

Site, Roller Clamp, Occlusion Clamp ● 84" Tubing

50 pcs/case

Priming Volume 57ml/Approx Inside Diameter 3.0 ± 0.05mm Outside Diameter 4.2 ± 0.05mm

SAFETY HUBER INFUSION SETS NON DEHP

Patient Safety & Comfort

- Minimizes damage to port septum
- Significantly reduces penetration and access forces
- Designed for patient comfort
- Easy-to-dress
- Helps maintain dressing integrity
- Prevents loss of drug potency due to DEHP

This device does not contain natural rubber latex.

Easy to Dress


- Ultra Low Profile
- Small Footprint
- Angled Tuning
- **Easy to Use**

Safety Features

- Built-in Safety Mechanism
- Smooth, Light Lockout
- Ergonomic, Integral Wing Design
- Audible "Click", Tactile Feel and visual indication of Lockout
- Needle Forward Design

Benefits

- Complies with OSHA standards
- Reduces the risk of needle stick injury
- Controlled, easy safety mechanism engagement
- Controlled, secure access
- Confirms safety mechanism activation
- Facilitates dual lumen port access
- Easy to view needle and site during access



SECURE TOUCH+ SAFETY HUBER INFUSION SET

7" tube, 'Y' injection site (right angle)

Catalog No.	Description	Hub Colo
37840S	19G x 1"	Brown
37841S	19G x 3/4"	Brown
37844S	20G x 1"	Yellow
37845S	20G x 3/4"	Yellow
37846S	20G x 1 1/2"	Yellow
37848S	22G x 1"	Black
37849S	22G x 3/4"	Black
37850S	22G x 1 1/2"	Black
20/hox		

7" tube (rig	ht angle)	
Catalog No.	Description	Hub Color
37854S	20G x 1"	Yellow
37855S	20G x 3/4"	Yellow
37856S	20G x 1 1/2"	Yellow
37858S	22G x 1"	Black
37859S	22G x 3/4"	Black
20/box		

HUBER INFUSION SETS



- This device does not contain natural rubber Latex
- Non Coring
- 90 degrees 316L stainless steel needle
- 'Y' injection site available
- o 7" microbore tubing
- A rotating wing providing less stress on both the tissue and the port.
- Polypropylene hubs, packaging color-coded according to gauge for easy Identification
- Occlusion clamp

NON DEHP
ROTATING WING

NON DEHP HUBER (NON -CORING) INFUSION SETS

7" tube, rotating wing,
'Y' injection site (right angle)

NON DEHP HUBER (NON -CORING) INFUSION SETS

7" tube, rotating wing

Catalog No	. Description	Color Code
27940R	19G x 1"	Brown
27941R	19G x 3/4"	Brown
27944R	20G x 1"	Yellow
27945R	20G x 3/4"	Yellow
27946R	20G x 1 1/2"	Yellow
27948R	22G x 1"	■ Black
27949R	22G x 3/4"	Black
27950R	22G x 1 1/2"	■ Black
50/box		

Catalog No	. Description	Color Code
27954R	20G x 1"	Yellow
27955R	20G x 3/4"	Yellow
27958R	22G x 1 "	■ Black
27959R	22G x 3/4"	■ Black
50/box		

HUBER NEEDLES

- This device does not contain natural rubber latex
- NON-DEHP
- Non Coring
- 316L stainless steel construction
- Available in both straight and 90 degree angles
- Needle volume of +0.5 ml
 Polypropylene hubs, color-coded according to gauge for easy
 Identification



HUBER NEEDLES

(non-coring) right angle

Catalog	No. Description	Hub Colo
26901	20G x 3/4"	Yellow
26902	22G x 3/4"	Black
26904	22G x 1"	Black
26095	19G x 3/4"	Brown
26906	22G x 1 1/2"	Black
26907	19G x 1"	Brown
26908	20G x 1"	Yellow
26909	20G x 1 1/2"	Yellow
50/box		



HUBER NEEDLES

(non-coring) straight

Catalog No.	Description	Hub Color
26920	20G x 2"	Yellow
26921	22G x 1"	Black
26922	22G x 1 1/2"	Black
26923	20G x 1 1/2"	Yellow
26924	20G x 1"	Yellow
26925	20G x 1 1/4"	Yellow
50/box		

SAFELET IV CATHETERS





- Easy dispenser pack
- Functional design casing permitting one action opening and no rolling
- Color-coded cap allows for easier identification of catheter size
- Translucent catheter hub and flashback chamber allows for easy detection of blood flashback at vein insertion
- Radio-opaque striped catheter
- Precision finished FEP catheter assures stable flow and eliminates catheter tip kink during venipucture

- Cannula with additional back bevel
- Can be connected to syringe by removing filter cap to expose Luer taper end
- Use of hydrophobic membrane filter eliminates blood leakage
- Close and smooth contact between catheter tip and inner needle enables safe and smooth venipuncture
- Crimped smooth engagement between catheter and hub, without use of adhesive, helps eliminate blood residuals and helps reduce blood coagulation or hemolysis

SPECIFICATIONS

50/ box, 200/case

		Stainless	Cannula Stainless Can	nula	
Item No.	Size	Gauge Outside I	Diameter Total Length	h Hub Color	
		(mı	m) (mm)		
26726	14G x 1 1/4"	14G 1.6	51.0	Orange	
26729	14G x 2"	14G 1.6	67.0	Orange	
26730	16G x 1 1/4"	16G 1.:	2 50.0	Medium Grey	
26731	16G x 2"	16G 1.:	2 66.0	Medium Grey	
26735	18G x 1 1/4"	18G 0.	9 47.0	Green	
26736	18G x 1 1/2"	18G 0.	9 54.0	Green	
26738	18G x 2"	18G 0.	9 66.0	Green	222G EXEL
26741	20G x 1"	20G 0.7	71 41.0	Pink Pink	EXEL
26742	20G x 1 1/4"	20G 0.7	71 48.0	Pink	SAFELET CATHETER SETUP PURPOSE PLANTON
26743	20G x 2"	20G 0.7	71 66.0	Pink	<u>~ 7.5</u> €]50
26746	22G x 1"	22G 0.5	56 40.0	Blue	
26751	24G x 3/4"	24G 0.4	11 34.0	Yellow	

MULTI-SAMPLE NEEDLES MULTI-SAMPLE LUER ADAPTERS MULTI-SAMPLE LUER HOLDERS

INJECTION PLUG

MULTI-SAMPLE (BLOOD-DRAW) NEEDLES

Catalog No. Description **Hub Color** 26500 18G x 1" Pink 26501 20G x 1" Yellow 26502 20G x 1 1/2" Yellow 26503 21G x1" Green 26504 21G x 1 1/2" Green 22G x 1" 26505 Black 22G x 1 1/2" 26506 Black 100/box, 1000/case

- This product does not contain natural rubber latex.
- Sterile
- Non-toxic
- Multi-sample needles permit several samples to be taken with a single puncture
- Sharp & smooth edges make penetration painless for patients and provides easy connection to rubber stoppers



LUER LOCK HOLDER WITH PRE-ATACHED MULTIPLE SAMPLE ADAPTER

Catalog No.

Description

26532

5-7.5ML HOLDER

50/box,200/case

Available in bulk pack

Features:

- Does not contain natural rubber latex.
- Sterile
- Individually packed.
- Luer Lock connection.
- Blue color thread.

MULTI-SAMPLE LUER ADAPTERS

MULTI-SAMPLE

LUER HOLDERS

Catalog No. Description

5-7.5ml holder

100/bag, 1000/case (Non-Sterile)

Catalog No. Description

26537

26530

Adapter

100 pcs/box, 1000 pcs/case





INJECTION PLUG

Catalog No.

Description

26539

Injection plug with cap LL

50/box,200/case



DENTAL NEEDLES



DENTAL NEEDLES PLASTIC HUB

 Catalog No.
 Description
 Length

 26555
 30G Short
 (21mm)

 26562
 30G Long
 (32mm)

 26557
 27G Long
 (32mm)

 100/box,1,000/case



Hub Color

Blue Blue

Yellow



Convenient
RED DOT
marks the
bevel position

NEW
Hub design
for greater needle
flexibility

- New hub design lowers the chances of needle breakage by 500%!
- New RED DOT on the needle hub marks the needle bevel position, even when the needle point is not visible
- Years of exacting experience in the production of hypodermic needles has been applied to the production of Dental Needles, virtually painless, atraumatic, and perfectly sharp to give maximum patient comfort. Dental Needle production is computer controlled in all processes from manufacture to inspection. In addition, every needle passes an optical inspection
- The cannula is siliconized
- Needles are sterilized by Ethylene Oxide and their sterility is guaranteed by testing every production lot prior to release
- Dental Needle hubs and boxes are color coded by gauge for easier identification of suitable needle and gauge
- Each needle is individually packaged, sterilized, and ready for immediate use
- Self-tapping ribs placed on the outside of the hub, and inside threads make for ease of use
- Larger cannula lumen (.0055")

SPINAL NEEDLES

- This product does not contain natural rubber latex.
- Convenient blister packaging
- Siliconized for smooth comfortable penetration
- Transparent plastic hub



'QUINKE' BEVEL SPINAL NEEDLES

Description	Hub Color
18G (1.20mm) x 3 ½"	Pink
20G (1.20mm) x 2 3/4"	Yellow
20G (0.90mm) x 3 ½"	Yellow
22G (0.70mm) x 2 3/4"	Black
22G (0.70mm) x 3 ½"	Black
25G (0.50mm) x 3 ½"	Orange
	18G (1.20mm) x 3 ½" 20G (1.20mm) x 2 ¾" 20G (0.90mm) x 3 ½" 22G (0.70mm) x 2 ¾" 22G (0.70mm) x 3 ½"

50/ box, 100/case

SURGICAL BLADES DISPOSABLE SCALPELS





- This product does not contain natural rubber latex.
- Made from premium quality Swiss
 Carbon Steel and British Stainless Steel
- Computer controlled grinding, sharpening and polishing, ensures consistently smooth incision properties
- Sterilized by GAMMA Radiation 2.5m -RADS
- A Vapor Corrosion Inhibitor (VCI) paper lining, leakage testing, and an aluminum foil peal pack help insure quality and integrity when ready for use

- Each production lot is sampled and conforms to ISO 2859, inspection level (S4) AQL 1.0 & 2.5
- Blade rigidity is tested to eliminate blade breakage
- Blades conform to ISO: 7740 (1995) & BS: 2982 (1986)
- Blades are ultrasonically cleaned prior to packing

STAINLESS STEEL SURGICAL BLADES

Catalog No.	Size
29500	10
29502	11
29504	12
29508	15
29515	22

100/box,1,000/case

CARBON STEEL SURGICAL BLADES

Catalog No.	Size
29526	11
29532	15

100/box,1,000/case

STAINLESS STEEL DISPOSABLE SCALPELS

Catalog No.	Siz
29550	10
29552	11
29556	15
29559	22


10/box,100/case

PRESSURE BANDAGES

ZorBand PRESSURE BANDAGES



- Sterile
- Safe
- Easy to use
- Increased efficiency
- Individually packaged
- More convenient than traditional dressings
- ZorBand reduces direct pressure time by over 40%
- Absorbent cellulose pad expands when contacted by blood
- Swelling action of pad creates a pressure compress
- Firmly applies pressure to puncture
- ZorBand tightly binds to patients' skin to guard from leakage and maintain pressure



Catalog No. Description

26834 Large 26836 X-Large 100/box, 1,000/case

BULK NON-STERILE SYRINGES

1cc BULK NS SYRINGE

Catalog No. Description Qty. Per Case


* BN 26048 Luer Slip 1000 * BN 26049 TB Luer Lock 4000 BN 26050 TB Luer Lock 4000

* Low Dead Space Plunger This device does not contain natural rubber latex.



Catalog No. Description Qty. Per Case

* BN 26200 Luer Lock 1500 * BN 26201 Luer Slip 1500



5-6cc BULK NS SYRINGE

Catalog No. Description Qty. Per Case

BN 26230 Luer Lock 2000 BN 26231 Luer Slip 2000

20-25cc BULK NS SYRINGE

Catalog No. Description Qty. Per Case

BN 26280 Luer Lock 1000 BN 26281 Luer Slip 1000

10-12cc BULK NS SYRINGE

Catalog No. Description Qty. Per Case

BN 26265 Luer Lock 1500 BN 26266 Luer Slip 1500



30-35cc BULK NS SYRINGE

Catalog No. Description Qty. Per Case

BN 26290 Luer Lock 400 BN 26291 Luer Slip 400 BN 26292 Cath Tip 400

50-60cc BULK NS SYRINGE

Catalog No. Description Qty. Per Case

 BN 26300
 Luer Lock
 400

 BN 26301
 Luer Slip
 400

 BN 26302
 Cath Tip with cap
 400



EXELINT INTERNATIONAL, CO.

Corporate Headquarters

2500 Santa Fe Avenue

Redondo Beach

California 90278 USA

Phone: (310) 649-0707

Fax: (310) 649-1178

E-mail: info@exelmed.com

Sales & Marketing

1510 Dr. MLK Jr. Street North

Saint Petersburg,

Florida 33704 USA

Phone: (800) 940-3935

Fax: (800) 308-5048

FOR INFO PLEASE CALL

1-800-940-EXEL(3935)

info@exelint.com

www.exelint.com