




Relief Pak® ColdSpot™ reusable cold packs

Relief Pak® reusable cold packs

**soothe discomfort from injury
packs are segmented for consistent
fill and convenient contour folding**

R ■ an effective cold therapy for sprains and strains, muscle soreness, bruises, swelling and tension headache


- maintains therapeutic cold for up to 30 minutes
- place pack in a chilling unit or freezer
- pliable and soft, even below freezing
- use with cold pack cover (sold below) or hand towel between skin and pack
- black urethane packs made of heavier material than blue vinyl and can be chilled to lower temperatures (not below 0°F)
- 1-year warranty



blue vinyl Relief Pak® ColdSpot™ packs

packs only		each	bulk (12 each)
L	circular (10" diameter)	11-1007 15.75	11-1007-12 170.00
M	standard (11" x 14")	11-1000 15.75	11-1000-12 170.00
N	neck (6" x 23")	11-1001 16.75	11-1001-12 180.00
O	oversize (11" x 21")	11-1002 28.50	11-1002-12 310.00
P	half size (7" x 11")	11-1003 9.50	11-1003-12 102.50
Q	quarter size (5" x 7")	11-1004 7.75	11-1004-12 85.00
R	throat (3" x 11")	11-1005 8.00	11-1005-12 85.00

black urethane Relief Pak® ColdSpot™ packs



packs only		each	bulk (12 each)
S	circular (10" diameter)	11-1257 16.25	11-1257-12 177.50
T	standard (11" x 14")	11-1250 16.25	11-1250-12 177.50
U	neck (6" x 23")	11-1251 17.00	11-1251-12 185.00
V	oversize (11" x 21")	11-1252 28.50	11-1252-12 310.00
W	half size (7" x 11")	11-1253 10.00	11-1253-12 110.00

Relief Pak® cold packs chilling units



- **front loading chilling unit:** 4.6 cubic foot capacity with 2 full width door shelves and 2 full width wire shelves. Adjustable thermostat control. Defrost drain and manual defrost. 39"H x 27"L x 22"D; 137 lb

X		11-0501 front loading	500.00
		11-0506 front loading with 12 Relief Pak® standard (11-1000)	635.00

- **top loading chilling unit:** 5.0 cubic foot capacity with 1 removable basket. Adjustable thermostat control, defrost drain and manual defrost. 32"H" x 30"L x 24"D; 110 lb

Y		11-0500 top loading	390.00
		11-0505 top loading with 12 Relief Pak® standard (11-1000)	510.00

Relief Pak® cold pack covers

covers only		each	bulk (12 each)
	circular (10" diameter)	11-1009 12.50	11-1009-12 135.00
	standard (13" x 16")	11-1010 12.50	11-1010-12 135.00
	neck (9" x 25")	11-1011 15.00	11-1011-12 162.50
	oversize (13" x 23")	11-1012 17.50	11-1012-12 190.00
	half/quarter (9" x 13")	11-1013 10.00	11-1013-12 110.00

*covers not shown