

- use for tip, key and palmar tests
- Baseline[®] hydraulic pinch gauge uses a hydraulic system to assure measurement accuracy, repeatability and product reliability
- therapist can support the pinch gauge during testing, yielding a more accurate result for all pinch tests
- manufactured and serviced in the USA
- results consistent with published Baseline[®] and Jamar[®] studies
- maximum reading remains until reset
- readings in lb and kg
- lightweight unit includes protective carrying case
- CE certified

50 lb standard capacity



- reliable
- 50 lb (23 kg) capacity
- 2½" diameter dial
- attractive red finish

LITE[®]
1 year warranty

- our most popular style
- 50 lb (23 kg) capacity
- 2½" diameter dial
- blue finish

Standard
2 year warranty

- maximum warranty
- 50 lb (23 kg) capacity
- 2½" diameter dial
- black finish

HD[®]
5 year warranty

12-0226 50 lb LiTE[®]

200.00

12-0235 50 lb standard

250.00

12-0222 50 lb HD[®]

300.00




- extra large 3½" diameter dial
- 50 lb (23 kg) capacity
- blue finish

Hi-Res[®]
2 year warranty

12-0239 50 lb HiRes[®]

275.00

extended range 100 lb capacity



- extended range 100 lb (45 kg) capacity
- extra large 3½" diameter dial
- blue finish

Hi-Res[®]
ER[™]
2 year warranty

12-0228 100 lb ER[™] HiRes[®]

300.00

- extended range 100 lb (45 kg) capacity
- electronic zero calibration system
- max clear button
- toggle maximum and current reading
- toggle lb / kg
- 2½" LCD gauge
- 2 AAA batteries included
- blue finish

digital
2 year warranty

12-0237 100 lb digital

650.00