

## Baseline® 5-position hydraulic pinch grip



use without paddle to accommodate small hands (same pinch width as standard pinch gauges)

5-grip width can accommodate large hands

- easy to adjust paddle accommodates to any size hand - other pinch gauges only have one pinch width
- pinch width without paddle is the same width as standard pinch gauge (mechanical and hydraulic)
- 5 pinch positions permits MVE and MMVE protocols to be used with tip, palmar and key tests

12-0480	5-position clinical pinch gauge	250.00
12-0490	5-position feedback pinch gauge	250.00
12-0480P	additional paddle for 5-position pinch gauge	25.00



- adjustable handle fits all size hands for accurate and repeatable pinch measurements
- 5-position strength test protocols (MVE and MMVE) now can be used for pinch measurements

### clinical model

- gauge faces away from subject (toward practitioner) to block readings from subject

12-0160 set w/case 560.00

#### set contents

12-0240	grip (200lb)	325.00
12-0480	pinch (50lb)	250.00
12-1014	goniometer	30.00

### feedback model

- gauge faces toward subject to enhance compliance and feedback

12-0161 set w/case 560.00

#### set contents

12-0271	grip (200lb)	325.00
12-0490	pinch (50lb)	250.00
12-1014	goniometer	30.00




feedback model

## Baseline® MMT (Manual Muscle Test) kit

**use your existing hand dynamometer for manual muscle testing!**

### here's how

- MMT handle converts ANY hydraulic grip strength (hand) dynamometer to an MMT (Manual Muscle Test) dynamometer
- MMT handle comes with 3 push pad attachments: straight, curved and small circular



**B+C MMT handle with 3 push pads**  
12-0580 fits all hydraulic dynamometers 150.00

**D grip strength dynamometer with MMT handle and 3 pads**

12-0590	LiTE®	1 year warranty	375.00
12-0591	standard	2 year warranty	475.00
12-0592	HD®	5 year warranty	575.00

**E grip strength dynamometer with MMT handle, 3 pads and stabilizer grip**

12-0593	LiTE®	1 year warranty	475.00
12-0594	standard	2 year warranty	575.00
12-0595	HD®	5 year warranty	675.00

**F dual-handled stabilizer grip only**  
12-0583 fits all hydraulic dynamometers 100.00

evaluation