

CanDo® individual exercise/rest mats




Polyethylene (PE) cross-link mats have 1.7 lb density

Polyurethane (PU) mats are IFD100

Enviro Safe® mats are environmentally friendly

color options



- therapy mat is available with or without center-fold
- mat is covered with 18 oz tough, fire-retardant phthalate-free PVC vinyl; cover is sewn around foam
- cover is strong, mildew-resistant, anti-bacterial and wipes clean with a damp cloth
- mat is made with square side walls and is reinforced at all stress points for added durability
- PE, PU, and EnviroSafe® Foam options:
  - polyethylene (PE) cross-link foam - 1 3/8" thick PE mats have a hard feel; used for more strenuous exercises
  - polyurethane (PU) foam - 2" thick PU mats have a softer feel; used when cushioning is most important
  - EnviroSafe® Foam - 2" thick softer feel, 1 3/8" thick medium firm
- choose color from chart on right; **royal blue** is default

center-folding mat

W x L (ft)	1 3/8" PE foam	2" PU foam	2" EnviroSafe®
2' x 4'	38-0500 50.00	38-0400 77.50	38-2400 77.50
2' x 5'	38-0501 60.00	38-0401 90.00	38-2401 90.00
2' x 6'	38-0502 65.00	38-0402 102.50	38-2402 102.50


non-folding mat

W x L (ft)	1 3/8" PE foam	2" PU foam	2" EnviroSafe®
2' x 4'	38-0510 45.00	38-0410 72.50	38-2410 72.50
2' x 5'	38-0511 52.50	38-0411 85.00	38-2411 85.00
2' x 6'	38-0512 57.50	38-0412 97.50	38-2412 97.50

EnviroSafe® foam mats are the eco-friendly alternative featuring:

- 100% toxin free, zero VOC's or off gas, hypoallergenic, latex-free, mold and mildew proof, 100% recyclable with no leaching chemicals. Conforms to California Fire Retardant Code 117.

CanDo® center-fold and non-folding mat



A non-folding mat with handles

W x L (ft)	1 3/8" PE foam	2" PU foam	1 3/8" EnviroSafe®	2" EnviroSafe®
4' x 4'	38-0315 150.00	38-0305 167.50	38-2315 152.00	38-2305 167.50
4' x 6'	38-0310 210.00	38-0300 225.00	38-2310 210.00	38-2300 225.00
4' x 7'	38-0316 235.00	38-0306 255.00	38-2316 235.00	38-2306 255.00
4' x 8'	38-0311 267.50	38-0301 282.50	38-2311 267.50	38-2301 282.50
5' x 7'	38-0312 272.50	38-0302 285.00	38-2312 272.50	38-2302 287.50
5' x 10'	38-0313 377.50	38-0303 392.50	38-2313 377.50	38-2303 392.50
6' x 8'	38-0317 377.50	38-0307 392.50	38-2317 377.50	38-2307 392.50
6' x 12'	38-0314 530.00	38-0304 565.00	38-2314 530.00	38-2304 565.00

B center-folding mat with handles

W x L (ft)	1 3/8" PE foam	2" PU foam	1 3/8" EnviroSafe®	2" EnviroSafe®
4' x 4'	38-0215 172.50	38-0205 182.50	38-2215 172.50	38-2205 182.50
4' x 6'	38-0210 230.00	38-0200 240.00	38-2210 230.00	38-2200 240.00
4' x 7'	38-0216 262.50	38-0206 272.50	38-2216 262.50	38-2206 272.50
4' x 8'	38-0211 287.50	38-0201 302.50	38-2211 287.50	38-2201 302.50
5' x 7'	38-0212 302.50	38-0202 325.00	38-2212 302.50	38-2202 325.00
5' x 10'	38-0213 407.50	38-0203 430.00	38-2213 407.50	38-2203 430.00
6' x 8'	38-0217 407.50	38-0207 430.00	38-2217 407.50	38-2207 430.00
6' x 12'	38-0214 557.50	38-0204 612.50	38-2214 557.50	38-2204 612.50

- reinforced handles can be used to carry and store mats
- see information about foam options (above)
- choose color from chart (above); **royal blue** is default

EnviroSafe® foam mats are the eco-friendly alternative with the following features:

- 100% toxin free, zero VOC's or off gas, hypoallergenic, latex-free, mold and mildew proof, 100% recyclable with no leaching chemicals. Conforms to California Fire Retardant Code 117.