

# FHC



**NEW**

**FHC7200**

## 5th WHEEL, 750 lb CAPACITY TRANSPORT STRETCHER

*Come see what our stretchers have to offer your facility. Future Health Concept's Transport Stretchers provide the best blend of features, value, and customer support in the industry.*


**CALL OR EMAIL US TODAY!**

[www.fhcusa.com](http://www.fhcusa.com)

Phone: (888)282-8644.

[sales@fhcusa.com](mailto:sales@fhcusa.com)


**RETRACTABLE CENTER FIFTH  
WHEEL PROVIDES STRAIGHT  
LINE MOBILITY ALONG WITH  
EFFORTLESS PIVOTING  
AND CORNERING**



**OMNI-DIRECTIONAL STEERING  
MAKES FOR EASY STORAGE  
IN TIGHT SPACES**




MADE IN THE USA




**TREND/REVERSE TREND +/-18°  
SIX BRAKE/STEER PEDAL LOCATIONS**




**DUAL FOWLER ASSISTED  
LIFT BACK REST WITH 90° OF  
ADJUSTABILITY**






**LARGE STORAGE AREA FOR  
PATIENT BELONGINGS**



**SIX I.V. POLE  
CONFIGURATIONS**



**ERGONOMICALLY  
LOCATED O<sup>2</sup>  
TANK HOLDER**



**Standard Features:**

- Retractable Fifth wheel provides effortless directional movement and cornering
- 750 lb weight capacity
- O2 storage under bed platform for better ergonomics. No bending over
- Easily store and organize fleets in tight spaces with 4 swiveling wheels
- Dual Cylinder Fowler lift assist on raised back rest
- 8" Tente total locking casters
- 6 Brake/steer pedals
- Dual pedestal lift mechanism
- Trendelenburg/reverse Trendelenburg +/-18°
- Lift/Trendelenburg pedals on both sides
- Fold-down side rails with built in push handles
- 4" deluxe mattress
- 2 hook removable telescopic I.V. pole
- 6 I.V. pole sockets located on all corners
- Large storage area underneath for patient items
- Full perimeter bumper to protect your walls and fixtures
- 4 restraint belt locations

**Optional Features:**

- Fold down push handles
- Manual knee gatch
- Manual Fowler