

FORA®

IR42 Non-Contact Forehead Thermometer




Key Features



High Accuracy with Advanced Infrared Non-Contact Technology

Measurement distance from 1.97" to 2.75" (5 to 7 cm)



3 Modes: Adults, Children (3-36 months), Surface



3 modes



Color Backlight



Modes	Lo (No Light)	Green Backlight	Red Backlight	Hi (No Light)
Adults	$\leq 89.5^{\circ}\text{F}$	89.6°F to 100.3°F	100.4°F to 109.4°F	$\geq 109.5^{\circ}\text{F}$
Children	$\leq 89.5^{\circ}\text{F}$	89.6°F to 99.8°F	99.9°F to 109.4°F	$\geq 109.5^{\circ}\text{F}$
Surface	$< 32^{\circ}\text{F}$	32°F to 212.2°F	NA	$> 212.2^{\circ}\text{F}$

Ergonomically Designed



One touch button

Benefits



- ✓ Get the reliable reading instantly
- ✓ Ergonomically Designed: Easy to use for everyone
- ✓ Check adults, children and surface temperature with one meter
- ✓ 30 Memories: Easy to trace temperature history
- ✓ Avoid cross infection

Specification





Item	Specification
Temperature Mode	Adult/Child/Surface
Temperature Unit	°F (Default) or °C
Forehead Measurement Range	89.6°F to 109.4°F (32°C to 43°C)
Surface Measurement Range	0°C to 100°C (32°F to 212°F)
Measurement distance	1.18" - 2.75" (3 to 7 cm)
Memory Capacity	30 measurements
Display Resolution	0.1°F/0.1°C
Acoustics	Normal temperature range: two short beeps Fever: one long and three short beeps
Backlight	Normal temperature range: Green Fever: Red
Dimension	6.12" (L) x 1.58" (W) x 1.55" (H)
Weight	0.14 lb (without battery)
Bluetooth	No
Power source	2 x AAA included