

# Noninvasive Blood Pressure Air Hoses


Comprehensive solutions for noninvasive blood pressure

## GE Air Hoses

GE Healthcare provides a comprehensive portfolio of air hoses for use with GE monitors connecting to GE CRITIKON® blood pressure (BP) cuffs. GE air hoses are offered in a variety of sizes and styles – including adult, pediatric and neonatal. GE air hoses are not made with natural rubber latex and help prevent air leakage by providing a secure connection to the BP cuff.



**Adult-Pediatric, Gray Shroud to Subminiature Connector**  
Compatibility: GE DINAMAP® Pro Series, GE DINAMAP ProCare, CARESCAPE® V100, Compact, MPS Series Monitors  
**Length: 3.6 m/12 ft    Part # 107363**



**Adult-Pediatric, Gray Shroud to Screw Connector**  
Compatibility: GE DINAMAP Pro Series, GE DINAMAP ProCare, CARESCAPE V100, Compact, MPS Series Monitors  
**Length: 3.6 m/12 ft    Part # 107365**  
**Length: 3.6 m/12 ft    Part # 107366**



**Adult-Pediatric, Gray Shroud to DINAClick® Connector**  
Compatibility: GE DINAMAP Pro Series, GE DINAMAP ProCare, CARESCAPE V100, Compact, MPS Series Monitors  
**Length: 3.6 m/12 ft    Part # 2058205-001**  
**Length: 7.2 m/24 ft    Part # 2058205-002**






**Adult-Pediatric, Round Connector to Subminiature Connector**

Compatibility: Eagle 4000, Tram® 100/200/250/300/350, Corometrics® 120


Length: 3.6 m/12 ft Part # 9461-203



**Neonatal, Round Connector to Female Slip Luer**

Compatibility: Eagle 4000, Tram 100/200/250/300/350


Length: 2.4 m/8 ft Part # 9461-205



**Adult-Pediatric, Round Connector to Screw Connector**

Compatibility: Eagle 4000, Tram 100/200/250/300/350, Corometrics 118

Length: 3.6 m/12 ft Part # 9461-214



**Neonatal, Round Connector to Female Slip Luer**

Compatibility: Eagle 4000, Tram 100/200/250/300/350

Length: 3.6 m/12 ft Part # 9461-215



**Adult-Pediatric, Screw Connector to Screw Connector**

Compatibility: GE DINAMAP 1846, 8100, 9300 and Plus Series Monitors

Length: 3.6 m/12 ft Part # 008841  
Length: 7.2 m/24 ft Part # 008842



**Adult-Pediatric, Screw Connector to Subminiature Connector**

Compatibility: GE DINAMAP 1846, 8100, 9300 and Plus Series Monitors

Length: 3.6 m/12 ft Part # 088845



**Neonatal, Gray Shroud to Neo-Snap\* Connector**

Compatibility: GE DINAMAP Pro Series, GE DINAMAP ProCare, CARESCAPE V100, Compact, MPS Series Monitors

Length: 3.6 m/12 ft Part # 2058206-001



**Neonatal, Gray Shroud to Female Slip Luer**

Compatibility: GE DINAMAP Pro Series, GE DINAMAP ProCare, CARESCAPE V100, Compact, MPS Series Monitors

Length: 3.6 m/12 ft Part # 107368



**Air Hose Instrument Connector Assembly**

Part # 330073

**For More Information**

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)

GE Healthcare  
8200 West Tower Avenue  
Milwaukee, WI 53223  
U.S.A.

[www.gehealthcare.com](http://www.gehealthcare.com)



© 2013 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

\*Trademarks of General Electric Company.

GE Healthcare, a division of General Electric Company.

GE Medical Systems *Information Technologies*, Inc., a General Electric company, doing business as GE Healthcare.