

Choose the Height and Width - Choose Your Drawer Configuration Choose Your Lock - Your Color - Your Accessories

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



MDS3030EKC06

Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

M-Series Features:

- Choice of lock options
- Full extension drawers slides for easy access
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant hard baked powder coat finish
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable



Engineered &
Manufactured
in the U.S.



MDS3027BC5

(upgraded to white anti-microbial breakaway bar and other options)

Base MDS30 Cart Dimensions:

- Inside Drawer Dims: 24.5" W x 17" D (62.2 x 43.2 cm)
- Cart Weight – 100-160 lbs, depending on model
- Shipping Dims: 56" H x 39" W x 26" D (142.2 cm x 99 cm x 66 cm)

The Choice is Yours!

X-Short
34" H x 36.75" W x 22" D
 32.875" Height of Work Surface
**21" of Configurable
 Drawer Space**



MDS3021E05

Short
37.25" H x 36.75" W x 22" D
 36.125" Height of Work Surface
**24" of Configurable
 Drawer Space**



MDS3024E05

Medium
40.5" H x 36.75" W x 22" D
 39.375" Height of Work Surface
**27" of Configurable
 Drawer Space**



MDS3027ELP06


Tall
43.75" H x 36.75" W x 22" D
 42.625" Height of Work Surface
30" of Configurable Drawer Space




MDS3030B06

Choose Your Drawer Configuration:


1. Choose a cart size (X-Short, Short, Medium and Tall)
2. Select **any combination** of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.




MDS3021



MDS3024



MDS3027



MDS3030

Choose Your Lock:

- **B** - Breakaway Lock
- **K** - Key-Lock
- **E** - Basic Electronic Pushbutton Lock
- **EKC** - Electronic Lock w/Keypad
- **ELP** - Electronic Lock w/Keypad and Proximity Reader
- **WKC** - Wireless Electronic Lock w/Keypad
- **WELP** - Wireless Electronic Lock w/Keypad and Proximity Reader
- View locking options online: <https://www.harloff.com/locking-options/>


Choose Your Color(s):

Available in any of 23 colors. An anti-microbial finish is available in Lemon Drop Yellow. View color options online: <https://www.harloff.com/finish-options/>

Choose Your Accessories:

Full line of drawer dividers and exterior accessories are available for all cart sizes. View accessory options online: <https://www.harloff.com/accessories/>

Three M-Series Widths to Choose From:





M-SERIES EMERGENCY PACKAGES

M-series

BASIC PACKAGE

MD30-EMG



Oxygen tank holder (O2HLDR)



Cardiac board and brackets (MD-CARDBRD)



IV pole (MD-IV2)



STANDARD PACKAGE

MD30-EMG1



Oxygen tank holder (O2HLDR)



Drop shelf (MD-DRPSHLF)



(2) Utility hooks (UHOOKS2)



(2 Sets) Standard 3" Drawer Dividers (MD30-DIV3-S)



Cardiac board and brackets (MD-CARDBRD)



IV pole (MD-IV2)



View similar products online at <https://www.harloff.com/product-category/m-series-accessories/>

Toll Free: (800) 433-4064 • www.harloff.com
medicalsales@harloff.com • Ph: (719) 637-0300

SUPER STAT PACKAGE

MD30-EMG3



Oxygen tank holder (O2HLDR)



Large Utility Basket (BASKETDM)



(2) Utility hooks (UHOOKS2)



Toprail, matching cart color (MD30-TOPRAIL)



Cardiac board and brackets (MD-CARDBRD)



Chart Holder (CHARTHLDRDM)



(1 Set) Standard 3" Drawer Dividers (MD30-DIV3-S)



(1 Set) Standard 6" Drawer Dividers (MD30-DIV6-S)



(1 Set) Standard 3" Drawer Tray Dividers (MD30-TRAYDIV3-S)




IV pole (MD-IV2)



Defibrillator Shelf (MD-DEFIB)



Locking Sharps Container (5qt) and Glove Box Holder (684811)



Suction shelf (SUCTIONSHLF)



Engineered & Manufactured in the U.S.



Color shown is Red. More than 20 additional colors are available.



View similar products online at <https://www.harloff.com/product-category/m-series-accessories/>

Toll Free: (800) 433-4064 • www.harloff.com
 medicalsales@harloff.com • Ph: (719) 637-0300