

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 \* Ph: (719)637-0300 Fax:(719)597-8273 \* www.harloff.com

# OptimAL Line Multi-Dose Bin Medication Cart, Standard Package #AL48MDBIN



#### **Construction:**

· Painted lightweight aluminum cabinet and drawers

### **Standard Accessories:**

- Internal narcotics box with Illinois Lock
- BEST® brand lock on cabinet
- Painted pull-out side shelf with recessed work area
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with durable work surface
- Ergonomic top rail/push handle
- Easy clean protective pontoon bumpers
- 5" ball-bearing premium rubber casters, all full swivel, 1 brake and 1 directional
- Hard-baked powder coat finish (Customer chooses color for cabinet)

## **Drawer Configuration:**

- 4 ea 4" drawers
- 1 ea 11" multipurpose storage drawer with locking narcotics box

## Thirteen available color choices:

Beige, Burgundy, Light Blue, Red, Teal, Yellow, Dark Blue, Light Gray, Mauve, Hammertone Grey Hammertone Blue, Hammertone Green and Sand

Colors may be viewed at www.harloff.com

## **Assembled Cart Dimensions:**

- 42" H x 31" W x 24" D (includes top rail)
- Work surface: 38" H (96.5 cm)