

Choose the Height and Width - Choose Your Drawer Configuration Choose Your Lock - Your Color - Your Accessories

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



MDS2424ELP05



Engineered &
Manufactured
in the U.S.



MDS2430B06

Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

M-Series Features:

- Choice of lock options
- Full extension drawers slides for easy access
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant hard baked powder coat finish
- Ergonomic side push handles can be mounted on either side
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable

Base MDS24 Cart Dimensions:

- Inside Drawer Dims: 18.5" W x 17" D (47 x 43.2 cm)
- Cart Weight – 95-135 lbs, depending on model
- Shipping Dims: 56" H x 33" W x 26" D (142.2 cm x 99 cm x 66 cm)

The Choice is Yours!

X-Short
34" H x 30" W x 22" D
 32.875" Height of Work Surface
**21" of Configurable
 Drawer Space**



MDS2421K14

Short
37.25" H x 30" W x 22" D
 36.125" Height of Work Surface
**24" of Configurable
 Drawer Space**



MDS2424E05

Medium
40.5" H x 30" W x 22" D
 39.375" Height of Work Surface
**27" of Configurable
 Drawer Space**



MDS2427EKC05

Tall
43.75" H x 30" W x 22" D
 42.625" Height of Work Surface
**30" of Configurable
 Drawer Space**



MDS2430KC5-18EH

Choose Your Drawer Configuration:


1. Choose a cart size (X-Short, Short, Medium and Tall)
2. Select **any combination** of 3", 6", 9" or 12" drawers that fit within your cart's configurable drawer space.

7 drawers x
 3" height = 21"




MDS2421

2 drawers x
 12" height = 24"



MDS2424

3 drawers x
 9" height = 27"



MDS2427

5 drawers x
 6" height = 30"



MDS2430

Choose Your Lock:

- **B** - Breakaway Lock
- **K** - Key-Lock
- **E** - Basic Electronic Pushbutton Lock
- **EKC** - Electronic Lock w/Keypad
- **ELP** - Electronic Lock w/Keypad and Proximity Reader
- **WKC** - Wireless Electronic Lock w/Keypad
- **WELP** - Wireless Electronic Lock w/Keypad and Proximity Reader
- View locking options online: <https://www.harloff.com/locking-options/>

Choose Your Color(s):

Available in any of 23 colors. An anti-microbial finish is available in Lemon Drop Yellow. View color options online: <https://www.harloff.com/finish-options/>

Choose Your Accessories:

Full line of drawer dividers and exterior accessories are available for all cart sizes. View accessory options online: <https://www.harloff.com/accessories/>

Three M-Series Widths to Choose From:

