

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 * Ph: (719)637-0300 Fax:(719)597-8273 * www.harloff.com

V-Series, Four Drawer Procedure Cart with Key Lock #V30-4K

Harloff's V-Series Carts are value-priced for use in a variety of hospital departments. This drawer configuration is often a popular choice for treatment, procedure and isolation applications.



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

https://www.harloff.com/finish-options/

Specifications

Construction:

Steel construction

Features:

- Key lock with two keys
- Aluminum side rail / push handle
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- V-Series accessories are customer installed
- Latex free
- Durable, hard-baked powder coat finish

Internal Drawer Dimensions:

• 23" W x 16.75" D (58.4 cm W x 42.5 cm D)!

Drawer Configuration:

- (1) 3" (15.2 cm)
- (3) 9" (22.9 cm)

Integrated Pull Handle on all drawers

Assembled Cart Dimensions

- 40" H x 29.5" W x 22" D (101.6 cm H x 74.9 cm W x 55.9 cm D)
- Writing surface: 40" (101.6cm)
- Cart weight: 111 pounds (50.5 kg)
- Approx. shipping weight: 155 pounds (70.5 kg)
- Approx. shipping dims: 57" h x 33" W x 26" D (144.8cm H x 83.82cm W x 66cm D)