



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Value Line 6 Drawer Cart with Key Lock, Standard Package # VLT6K

Harloff's value medication carts offer many of the same features our customers love in our large med carts at an economical price point. This economy med cart has become our fastest growing small med cart.

Specifications



Construction:

- Painted steel cabinet and drawers

Features:

- Tall cabinet with key locking system with 4 keys
- Painted pull-out side shelf
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing premium casters, all full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Hard-baked powder coat finish
- Certified "Latex Free"
- Twelve year manufacturer's warranty

Drawer Configuration:

- 4 - ea 3" multipurpose storage drawers
- 2 - ea 9.75" drawers

Internal Drawer Dimensions:

- 15" W x 16.75" D (38.1 cm x 42.5 cm)

Assembled Cart Dimensions:

- 47" H x 21.3" W x 21.5" D
(119.4 cm H x 54.1 cm W x 54.6 cm D)
- Writing surface: 43.25" (109.9 cm)
- Cart net weight: 110 pounds (50 kg)
- Approx. shipping weight: 155 pounds (70.45 kg)

Thirteen available color choices:

Beige, Burgundy, Light Blue, Dark Blue, Light Gray, Mauve, Red, Teal, Yellow, Hammertone Blue, Hammertone Gray, Hammertone Green and Sand