

Common Features

PARALLEL BARS – MODELS 134X,
1386, 1387, 1388, 1389:



- 1½" diameter stainless steel handrails.
- Heavy gauge steel base plates, uprights and fittings. Black powder coated protective finish.
- Each upright telescopes in 1½" increments and locks into (10) height positions with fail-safe ergonomic locking pin.
- Width adjusts with triangular control knobs on each upright.

A. Model 1386, 1387 Bariatric Parallel Bars, Height & Width Adjustable

SPECIAL FEATURES:



- Designed to accommodate Bariatric patients. 600 lb. weight capacity.
- 10 ft. long gray laminate platform with anti-slip surface treads and tapered ends for wheelchair use.
- (6) steel uprights.
- Width adjusts from 22" to 35".
- 38" standard clearance between uprights.

MODEL	LENGTH	WIDTH*	HEIGHT
1386	10'	22" - 35"	29" - 42"
1387 (not shown)	7'	22" - 35"	29" - 42"

*Platform is 45" wide.

B. Model 1388, 1389 Floor Mounted Bariatric Parallel Bars

SPECIAL FEATURES:



- (6) floor-mounted uprights with mounting plates. Floor fasteners not included.
- Width adjusts from 22" to 35".

MODEL	LENGTH	WIDTH	HEIGHT
1388	10'	22" - 35"	29" - 42"
1389 (not shown)	7'	22" - 35"	29" - 42"

C. Series 1340 Floor Mounted Parallel Bars



SPECIAL FEATURES:

- (4) Floor-mounted uprights with mounting plates. Floor fasteners not included.
- (6) uprights for bars 14' or longer.
- Width adjusts from 15" to 28" wide.

MODEL	LENGTH	WIDTH	HEIGHT
1340	10'	15" - 28"	29" - 42"
1341 (not shown)	7'	15" - 28"	29" - 42"
1342 (not shown)	12'	15" - 28"	29" - 42"
1343 (not shown)	14'	15" - 28"	29" - 42"
1344 (not shown)	16'	15" - 28"	29" - 42"
1345 (not shown)	18'	15" - 28"	29" - 42"
1346 (not shown)	20'	15" - 28"	29" - 42"

A. Model 1386



B. Model 1388



C. Model 1340

