


A. Model 1440-57-113



## Common Features

DUAL LIFT POWERMATIC® MAT PLATFORMS — SERIES 1440, 1449:



- Power height adjustment from 19" wheelchair access height to 28"H.
- Two direct drive, synchronized lift mechanisms for unsurpassed power.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" High-Density urethane foam padding for comfort and rounded corners for safety.
- "Hands-Free" foot control.
- Hospital grade plug.
- 110V, 60Hz, 2.8 Amps. U.L. Listed

## A. Series 1440 Dual Lift Powermatic® Mat Platforms




MODEL	WIDTH	LENGTH	HEIGHT
1440-46-113	4'	6'	19" - 28"
1440-46-119	4'	6'	19" - 28"
1440-47-113	4'	7'	19" - 28"
1440-47-119	4'	7'	19" - 28"
1440-57-113	5'	7'	19" - 28"
1440-57-119	5'	7'	19" - 28"
1440-68-113	6'	8'	19" - 28"
1440-68-119	6'	8'	19" - 28"

*You must specify upholstery color with order.*

### COLORS:

**Laminate** – Black laminate base with satin aluminum trim and levelers.

### UPHOLSTERY – 113 BLUE OR 119 GRAY



### SPECIAL ORDER TOP COLORS



Model 1440-47-119

