

Common Features

BLOOD DRAWING CHAIRS — MODELS 2191, 2196:

- Designed to save time and increase efficiency.
- Gray powder-coated steel frame.
- One-piece contoured Black polypropylene seat 16"W x 16"D x 18"H.
- Armrests are interchangeable and can be set up for left or right hand use.
- One Armrest pivots 360°. The other Armrest is non-pivoting.
- Both Armrests have 10" height range.
- Easy-to-clean Gray laminate finish.
- Shipped ready-to-assemble.

A. Model 2191 Blood Chair – Standard Height



SPECIAL FEATURES:

- Seat is comfortable 19" height.
- Both Armrests adjust from 26" to 36" high with control knobs.

MODEL	LENGTH	DEPTH	HEIGHT
2191	29"	18"	19"

B. Model 2196 Blood Chair with Drawer



SPECIAL FEATURES:

- Same as Model 2191 but equipped with (1) polystyrene supply drawer. 8"W x 19"D x 3½"H.

MODEL	LENGTH	DEPTH	HEIGHT
2196	37½"	18"	19"

A. Model 2191



B. Model 2196



C. Model 2198



C. Model 2198 Bariatric Blood Chair



SPECIAL FEATURES:

- Extra-wide 24½" x 18½" deep x 22" high Black vinyl padded seat and backrest.
- "Flip-up" L-shaped, padded right hand Armrest folds upright 100° for convenient patient access. Folds down flat in use.
- Left Armrest is non-folding. Armrests are not interchangeable. Only available as shown.
- Both armrests are adjustable in height from 27" to 36" high with easy-to-use control knobs.
- Gray powder coated steel legs and frame
- Black vinyl is stocking color.

MODEL	WIDTH	DEPTH	HEIGHT
2198	34"	29"	37½"

D. Models 2000, 2002

Foot Stools



- Durable and attractive Silver Vein powder-coated steel finish.
- Anti-slip treaded 14" x 11" step.
- Metal tubular legs with safety tips.

MODEL	LENGTH	WIDTH	HEIGHT
2000	14"	11"	9"
2002	Same but packed (2) per carton		

D. Model 2000



E. Models 2030, 2032

Foot Stools with Safety Handrail



- Same as Model 2000 but with handrail.
- 33"H steel handrail snaps in without tools.

MODEL	LENGTH	WIDTH	HEIGHT
2030*	14"	11"	33"
2032	Same but packed (2) per carton		

*Platform is 9" high

E. Model 2030



F. Models 2216-346, 2216-927 Laminate Foot Stools



- Heavy duty, laminate wood construction in 346 Oak or 927 Folkstone Gray.
- Large 20" x 16" non-slip tread step.

MODEL	LENGTH	WIDTH	HEIGHT
2216-346	20"	16"	8"
2216-927	20"	16"	8"

F. Model
2216-927



F. Model
2216-346



SEE BARIATRIC FOOT STOOLS ON PAGE 34.