

## Common Features

**BARIATRIC FOOT STOOLS —  
MODELS 2018, 2038**



- Heavy duty 1000 lbs. weight capacity.
- Durable, maintenance-free, powder-coated Silver tone finish.
- Footprint of anti-tip base is 25" x 20".

### A. Model 2018 Bariatric XL Safe Step

**SPECIAL FEATURES:**

- Light weight for easy use – ONLY 6½ lbs.
- All-aluminum welded construction with Anti-Tip wide base.
- Anti-Slip treaded, extra-large 19" x 15" x 8"H platform with rubber angled leg tips.
- MRI SAFE...non-ferrous and non-magnetic.

MODEL	LENGTH	WIDTH	HEIGHT
2018	19"	15"	8"

### A. Model 2018



### B. Model 2038 Bariatric XL Safe Step with Handrail

**SPECIAL FEATURES:**

- Weighs only 10½ lbs. for ease of use.
- Heavy duty tubular aluminum 15" x 19" x 8" high platform with unique anti-tip design.
- 36" High sturdy steel safety handrail snaps in without tools.
- Anti-slip gripper strips & rubber angled leg tips.

MODEL	LENGTH	WIDTH	HEIGHT
2038	19"	15"	36"*

\*Platform is 8" high

### B. Model 2038



### C. Model 2010



### C. Model 2010 Bariatric Footstool

- Durable 1" bent steel tubing in maintenance-free Silver Vein finish with safety tips.
- 9" High non-skid "anti-slip" ribbed platform.
- Legs are cross-braced to accommodate up to 500 lbs.

MODEL	LENGTH	WIDTH	HEIGHT
2010	14"	10"	9"

### D. Model 2218



### D. Model 2218 Bariatric 2-Step Laminate Footstool

- Heavy-duty 500 lb. capacity wood construction in Natural Oak laminate.
- (2) 4" high steps with safety "anti-slip" treads. First step is 12" deep; top step is 16" deep. Bottom step is detachable if necessary.
- BriteStep 3mm yellow edging.

MODEL	LENGTH	WIDTH	HEIGHT
2218	28"	24"	8"



### E. Model 4074

### E. Model 4074 Extra-Wide Bariatric Treatment Table



- 3" extra-thick High-Density urethane foam top for bariatric use.
- Extra-wide 34" x 78" long upholstered top.
- Reinforced Oak laminated wood construction can accommodate patients up to 800 lbs.
- (3) Storage cabinets, each with laminate interiors and (1) adjustable shelf.
- Doors have concealed, self-closing, Euro-hinges.

**COLORS:**

Laminate – Natural Oak  
Upholstery – 17 standard vinyl colors.  
(FOLD OUT FRONT COVER)

MODEL	LENGTH	WIDTH	HEIGHT
4074	78"	34"	30"

**OPTIONS FOR MODEL 4074**

27 – Paper Dispenser & Cutter Combo