

A. Model 2153

Economy Air-Lift Stool

- Single lever height Control Handle.

MODEL	DIAMETER	HEIGHT
2153	16"	20"- 27"
2153-S <i>New Product!</i>	16"	17"- 22"



B. Model 2157

Economy Air-Lift Stool

- 360° circular ring Control Handle.

MODEL	DIAMETER	HEIGHT
2157	16"	20"- 27"
2157-S <i>New Product!</i>	16"	17"- 22"



Common Features

"CHOICE LINE™" AIR-LIFT STOOLS MODELS 2153, 2157:

- Hand operated air-spring pneumatic height adjustment.
- 16" diameter seat, 3½" thick foam.
- (5) 2" Twin wheel black nylon casters.
- Durable 5-leg, black composite base.

COLORS: (FOLD OUT FRONT COVER)

Upholstery -17 standard vinyl colors.



A. Model 2153



B. Model 2157-S

C. Model 2151

Heavy-Duty Air-Lift Stool

- Heavy-duty Pneumatic Lift with "D" chrome Control Handle.
- 16" Diameter Seat, 3½" thick foam seat.
- 4" stroke cylinder.
- Durable Aluminum 5-Leg Base.
- (5) 2½" Twin wheel gray nylon casters.
- Heavy-duty 400 lb. weight capacity.
- Available in 17 standard vinyl colors.



MODEL	DIAMETER	HEIGHT
2151	16"	20"-24"

D. Model 2154

"Dual" Foot & Hand Operated Stool

- Pneumatic Stool has BOTH Foot & Hand Controls for foot or hand operation.
- 5" stroke cylinder.
- Includes 360° ring Hand Control & circular Metal Foot Control.
- Durable 5-leg Polished Aluminum Base .
- (5) 2" Twin wheel black nylon casters.
- Available in 17 standard vinyl colors.



MODEL	DIAMETER	HEIGHT
2154	16"	19"-24"

C. Model 2151



D. Model 2154



E. Model 2156-707

Economy Air-Lift Stool with Backrest

- Padded backrest, 10" high x 15" wide.
- Air-spring pneumatic height adjustment.
- 360° circular ring Control Handle.
- Durable 5-leg, black composite base.
- 16" diameter seat, 3½" thick foam.
- (5) 2" Twin wheel black nylon casters.
- Black vinyl upholstery.

MODEL	DIAMETER	HEIGHT
2156-707	16"	20"- 27"

E. Model 2156-707

