


**B. Model 1449-57**  
*Shown with option:*  
**3002-57 – Mat**



**C. 3002-57 Mat attached with**  
**VELCRO® brand fasteners**

**B. Series 1449**  
**Dual Lift Powermatic® Mat**  
**Platforms with Black Laminate Tops**



- Will accommodate 2" thick reversible, removable mat, ordered separately. See #3002 below.
- Black laminate rimless subtop.
- Power height adjustment from 19" wheelchair access height to 28"H.
- Two direct drive synchronized lift mechanisms for unsurpassed power.
- "Hands-free" foot control.
- Hospital grade plug.
- 110V, 60Hz, 2.8 Amps. U.L. Listed

MODEL	WIDTH	LENGTH	HEIGHT	OPT.PAD
1449-46	4'	6'	19"-28"	3002-46
1449-47	4'	7'	19"-28"	3002-47
1449-57	5'	7'	19"-28"	3002-57
1449-68	6'	8'	19"-28"	3002-68

**C. Series 3002**  
**Platform or Floor Mats**



- Nylon-reinforced vinyl mat with 2" thick 100 I.L.D. urethane foam padding.
- (4) handles for easy moving and hanging.

**COLORS:**

**Reversible** – Blue on one side, Gray on the other with a Gray 2" border.

MODEL	WIDTH	LENGTH	HEIGHT
3002-46	4'	6'	2"
3002-47	4'	7'	2"
3002-57	5'	7'	2"
3002-68	6'	8'	2"

**D. Series 3042**  
**Folding Floor Mats**




- Accordion folding, 2' wide panels.
- 2" extra-firm foam.
- Nylon-reinforced vinyl cover.
- Easily joined to other mats with Velcro® brand fasteners on mat ends.
- Two-tone Blue and Gray.

MODEL	WIDTH	LENGTH	HEIGHT
3042-46	4'	6'	2"
3042-48	4'	8'	2"



**C. Model 3002-47**



**D. Model 3042-46**

# MAT PLATFORMS

A. Model 1441-47-113



## A. Series 1441 Econo Electric Hi-Lo Mat Platforms



**FEATURES:**

- Power height adjustment from 20" wheelchair height to 29"H.
- Two direct drive, synchronized lift mechanisms for greater power. 600 lbs. patient lift capacity.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" High-Density urethane foam padding for comfort and rounded corners for safety.

- Sturdy black laminate base with leveling legs.
- "Hands-Free" foot control.
- 110V, 60Hz, 1.3 Amps. Hospital grade plug.

MODEL	WIDTH	LENGTH	HEIGHT
1441-47*	4'	7'	20"-29"
1441-57*	5'	7'	20"-29"

**UPHOLSTERY - 113 BLUE OR 119 GRAY**



**SPECIAL ORDER TOP COLORS**



B. Model 1425-68-119



## B. Series 1425 Crank Hydraulic Mat Platforms



**FEATURES:**

- Durable, super-strength crank hydraulic system can lift 750 lbs. from 20" wheelchair height to 30" high.
- Easy-to-use hand crank height adjustment with fold away handle.
- Hydraulic cylinders are built into each (4) corner steel leg.
- Seamless reinforced vinyl is tear resistant.
- Steel frame is durable black powder coated finish.

- 2" High-Density urethane foam top.
- Legs have levelers for uneven floors.
- 4-Leg "Open" design is patient sling-lift accessible on all four sides.

MODEL	WIDTH	LENGTH	HEIGHT
1425-46*	4'	6'	20"-30"
1425-47*	4'	7'	20"-30"
1425-57*	5'	7'	20"-30"
1425-68*	6'	8'	20"-30"

*Shown with optional Bolsters (see page 46)*

*\*You must specify upholstery color with order.*