



# QUALITY-LINE TREATMENT TABLES



H-Brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.

HEIGHT CLEARANCE UNDER H-BRACE IS 7" WITH 31" HIGH TABLE

MODEL 4002



MODEL 4024



MODEL 4041  
Shown with option: 05 - Headrest

## A. Series 4002 H-Brace Table

- 4002\*: 72"L x 24"W x 31"H
- 4002-027: 72"L x 27"W x 31"H
- 4002-030\*: 72"L x 30"W x 31"H
- 4002-030-078: 78"L x 30"W x 31"H

## B. Series 4024 H-Brace Table with Laminate Storage Shelf.

- 4024\*: 72"L x 24"W x 31"H
- 4024-027: 72"L x 27"W x 31"H
- 4024-030\*: 72"L x 30"W x 31"H
- 4024-030-078: 78"L x 30"W x 31"H

## C. Series 4041 H-Brace Table with Laminate Shelf and Open-Shelf Storage Cabinet.

- 4041: 72"L x 24"W x 31"H
- 4041-027: 72"L x 27"W x 31"H
- 4041-030: 72"L x 30"W x 31"H
- 4041-030-078: 78"L x 30"W x 31"H

## D. Series 4043 Table with All-Laminate Storage Cabinet with Sliding Doors.

- Interior has divider and (1) adjustable shelf.
- 4043: 72"L x 24"W x 31"H
- 4043-027: 72"L x 27"W x 31"H
- 4043-030\*: 72"L x 30"W x 31"H
- 4043-030-078: 78"L x 30"W x 31"H



MODEL 4043

\*Same table available available as a Green-Line Table.

## Common Features

QUALITY LINE TREATMENT TABLES — SERIES 4002, 4024, 4041, 4043:

- 2" High-Density urethane foam top.
- Resist-All™ vinyl upholstery finish.

## COLORS:

Laminate – Natural Oak, Wild Cherry, Folkstone Gray  
Upholstery – 18 standard vinyl colors

400 LBS

## OPTIONS

FOR 4000 SERIES TABLES

- \*00 – Special Table Height (N/A Model 4043)
- \*04 – Table Safety Strap
- \*05 – Headrest, adjustable
- \*06 – Backrest, adjustable
- \*07 – Nose Cutout
- \*08 – Drawer (N/A with option 69)
- \*09 – Armboard, adjustable
- \*15 – Safety Rail (only Model 4043)
- \*16 – Stirrups, adjustable
- \*19 – Slide Shelf, laminate
- 23 – Paper Dispenser
- 25 – Paper Cutter
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small
- \*69 – Gas Spring Backrest
- \*99L – Left Hand Table Model

350 LBS

\*FACTORY INSTALLED

## ENHANCED LEG HARDWARE SYSTEM



Attach or detach legs quickly and easily with heavy-duty 5/16" bolts into embedded metal threaded inserts for durable connection. Easily tighten and re-tighten as needed.