

**A. Model 4795**  
**Bariatric Hi-Lo Treatment**  
**Table with Power Backrest**



**FEATURES:**

- Extra-wide 32" top accommodates larger patients up to 600 Lbs.
- Power height motorized elevation from 21" wheelchair height to 30" high.
- Two direct-drive synchronized lift actuators provide unsurpassed power.
- Motorized Power Backrest elevates 0° to 75° for ease of patient positioning.
- Convenient Hand Control operation.
- 2" High-Density urethane foam top for patient comfort.
- Black laminate frame and base plate, with aluminum trim and levelers.
- Hospital grade plug.
- 110V, 60Hz, 3.5 Amps. U.L. Listed.

**COLORS:**

**Upholstery** – 17 standard vinyl colors.  
*(FOLD OUT FRONT COVER)*

MODEL	LENGTH	WIDTH	HEIGHT
4795	78"	32"	21" - 30"

**OPTIONS**

FOR MODELS 4086, 4795

27 – Paper Dispenser & Cutter Combo

SEE PAGE 4



**A. Model 4795**



**Model 4795 shown in flat position**

**B. Model 4086**  
**Bariatric Treatment Table**



**FEATURES:**

- (6) high pressure laminate legs and a sturdy plywood sub-top.
- H-Brace with full width laminate shelf. 4" clearance under shelf.
- 2" High-Density urethane foam top.

**COLORS:**

**Laminate** – Natural Oak

**Upholstery** – 17 standard vinyl colors.  
*(FOLD OUT FRONT COVER)*

MODEL	LENGTH	WIDTH	HEIGHT
4086	78"	30"	26"

For 31" Height model – see A9087 (page 87)



**B. Model 4086**