



A. Model 4750

New Design!

Common Features

POWERMATIC® TREATMENT TABLES — MODELS 4750, 4756, 4757:



- 2" High-Density urethane foam vinyl top.
- "Hands Free" foot control.
- Includes Paper Dispenser and Cutter.
- Hospital grade plug.
- 110V, 60Hz, 3.5 Amps. U.L. Listed.
- New contoured black base plate with (4) levelers and satin aluminum trim.

COLORS:

Laminate – Folkstone Gray lift enclosure.

Upholstery – 17 standard vinyl colors.

(FOLD OUT FRONT COVER)

A. Model 4750 Powermatic® Table with Flat Top



SPECIAL FEATURES:

- Power height adjustment from 19" wheelchair height to 31" high.

MODEL	LENGTH	WIDTH	HEIGHT
4750	76"	27"	19"- 31"

B. Model 4756 Powermatic® Table with Gas-Spring Backrest

SPECIAL FEATURES:

- Power height adjustment from 19" wheelchair height to 31".
- Gas-spring pneumatic 30" backrest raises to 75° with (1) lever control.

MODEL	LENGTH	WIDTH	HEIGHT
4756	76"	27"	19"- 31"

Model 4757 (not shown) Powermatic® Table with Power Backrest

- Same as Model 4756 except Power 30" backrest elevates to 75° with foot control.

MODEL	LENGTH	WIDTH	HEIGHT
4757	76"	27"	19"- 31"



B. Model 4756

New Design!

OPTIONS

FOR MODELS 4750, 4756, 4757, 4766, 4767

- *04 – Table Safety Strap
- *07 – Nose Cutout (4750 only)
- 32 – Pillow, full size
- 35 – Pillow, small
- *99L – Left Hand Table Model

*FACTORY INSTALLED

SEE PAGE 4